

CX

Generated by Doxygen 1.8.6

Thu Feb 20 2014 01:26:26

Contents

1	Main Page	1
1.1	Installation	1
1.2	Examples and Tutorials	2
2	Modules	3
3	Blocking Code	3
4	license	3
5	Module Index	3
5.1	Modules	3
6	Namespace Index	4
6.1	Namespace List	4
7	Hierarchical Index	4
7.1	Class Hierarchy	4
8	Class Index	5
8.1	Class List	5
9	File Index	7
9.1	File List	7
10	Module Documentation	8
10.1	Data	8
10.1.1	Detailed Description	8
10.2	Entry Point	9
10.2.1	Detailed Description	9
10.2.2	Function Documentation	9
10.2.3	Variable Documentation	9
10.3	Error Logging	10
10.3.1	Detailed Description	10
10.3.2	Enumeration Type Documentation	10
10.4	Input Devices	11
10.4.1	Detailed Description	11
10.5	Randomization	12
10.5.1	Detailed Description	12

10.6 Sound	13
10.6.1 Detailed Description	13
10.7 Timing	14
10.7.1 Detailed Description	14
10.7.2 Variable Documentation	14
10.8 Utility	15
10.8.1 Detailed Description	15
10.9 Video	16
10.9.1 Detailed Description	16
11 Namespace Documentation	17
11.1 CX Namespace Reference	17
11.1.1 Detailed Description	18
11.1.2 Enumeration Type Documentation	18
11.2 CX::Draw Namespace Reference	18
11.2.1 Detailed Description	19
11.2.2 Function Documentation	19
11.3 CX::Instances Namespace Reference	20
11.3.1 Detailed Description	20
11.4 CX::Private Namespace Reference	20
11.4.1 Detailed Description	21
11.5 CX::Util Namespace Reference	21
11.5.1 Detailed Description	22
11.5.2 Function Documentation	22
12 Class Documentation	26
12.1 CX::CX_Clock Class Reference	26
12.1.1 Detailed Description	26
12.1.2 Member Function Documentation	26
12.2 CX::Util::CX_CoordinateConverter Class Reference	27
12.2.1 Detailed Description	28
12.2.2 Constructor & Destructor Documentation	28
12.2.3 Member Function Documentation	28
12.3 CX::CX_DataFrame Class Reference	30
12.3.1 Detailed Description	31
12.3.2 Member Function Documentation	31
12.4 CX::CX_DataFrameCell Class Reference	37
12.4.1 Detailed Description	38

12.4.2	Constructor & Destructor Documentation	39
12.4.3	Member Function Documentation	39
12.5	CX::CX_DataFrameColumn Class Reference	40
12.5.1	Detailed Description	41
12.6	CX::CX_DataFrameRow Class Reference	41
12.6.1	Detailed Description	41
12.7	CX::Util::CX_DegreeToPixelConverter Class Reference	41
12.7.1	Detailed Description	41
12.7.2	Constructor & Destructor Documentation	42
12.7.3	Member Function Documentation	43
12.8	CX::CX_Display Class Reference	43
12.8.1	Detailed Description	44
12.8.2	Member Function Documentation	44
12.9	CX::CX_FinalSlideFunctionInfo_t Struct Reference	47
12.9.1	Detailed Description	48
12.10	CX::CX_InputManager Class Reference	48
12.10.1	Detailed Description	48
12.10.2	Member Function Documentation	48
12.11	CX::CX_Joystick Class Reference	50
12.11.1	Detailed Description	50
12.11.2	Member Function Documentation	50
12.12	CX::CX_JoystickEvent_t Struct Reference	51
12.12.1	Detailed Description	52
12.12.2	Member Enumeration Documentation	52
12.13	CX::CX_Keyboard Class Reference	52
12.13.1	Detailed Description	53
12.13.2	Member Function Documentation	53
12.14	CX::CX_KeyEvent_t Struct Reference	53
12.14.1	Detailed Description	53
12.14.2	Member Enumeration Documentation	54
12.14.3	Member Data Documentation	54
12.15	CX::Util::CX_LengthToPixelConverter Class Reference	54
12.15.1	Detailed Description	54
12.15.2	Constructor & Destructor Documentation	54
12.15.3	Member Function Documentation	55
12.16	CX::CX_Logger Class Reference	55
12.16.1	Detailed Description	56

12.16.2 Member Function Documentation	56
12.17CX::CX_Mouse Class Reference	58
12.17.1 Detailed Description	58
12.17.2 Member Function Documentation	58
12.18CX::CX_MouseEvent_t Struct Reference	59
12.18.1 Detailed Description	59
12.18.2 Member Enumeration Documentation	60
12.19CX::CX_RandomNumberGenerator Class Reference	60
12.19.1 Detailed Description	61
12.19.2 Constructor & Destructor Documentation	61
12.19.3 Member Function Documentation	61
12.20CX::CX_SafeDataFrame Class Reference	67
12.20.1 Detailed Description	67
12.21CX::CX_Slide_t Struct Reference	68
12.21.1 Detailed Description	68
12.21.2 Member Enumeration Documentation	68
12.21.3 Member Data Documentation	68
12.22CX::CX_SlidePresenter Class Reference	69
12.22.1 Detailed Description	70
12.22.2 Member Function Documentation	70
12.23CX::CX_SlideTimingInfo_t Struct Reference	73
12.23.1 Detailed Description	73
12.23.2 Member Data Documentation	73
12.24CX::CX_SoundObject Class Reference	73
12.24.1 Detailed Description	74
12.24.2 Member Function Documentation	74
12.25CX::CX_SoundObjectPlayer Class Reference	77
12.25.1 Detailed Description	78
12.25.2 Member Function Documentation	78
12.26CX::CX_SoundStream Class Reference	79
12.26.1 Detailed Description	79
12.26.2 Member Function Documentation	79
12.27CX::CX_SoundStreamConfiguration_t Struct Reference	83
12.27.1 Detailed Description	83
12.27.2 Member Data Documentation	83
12.28CX::CX_SP_Configuration Struct Reference	84
12.28.1 Detailed Description	85

12.29CX::CX_SP_PresentationErrors_t Struct Reference	85
12.29.1 Detailed Description	85
12.30CX::Util::CX_TrialController Class Reference	85
12.30.1 Detailed Description	86
12.30.2 Member Function Documentation	86
13 File Documentation	87
13.1 advancedChangeDetectionTask.cpp File Reference	87
13.1.1 Detailed Description	88
13.2 basicChangeDetectionTask.cpp File Reference	88
13.2.1 Detailed Description	88
13.3 nBack.cpp File Reference	89
13.3.1 Detailed Description	89
Index	90

1 Main Page

ofxCX (hereafter referred to as [CX](#)) is a "total conversion mod" for openFrameworks (often abbreviated oF) that is designed to be used for creating psychology experiments.

The most well-organized way to access the documentation is go to the [Modules](#) page. The best way to get an overview of how [CX](#) works is to look at the [Examples and Tutorials](#).

To learn about presenting visual stimuli, go to the [Video](#) page.

To learn about auditory stimuli, go to the [Sound](#) page.

To learn how to store and output experiment data, see the [Data](#) page.

To learn about random number generation, see the [Randomization](#) page.

To learn about how [CX](#) logs errors and other runtime information, see the [Error Logging](#) page.

1.1 Installation

In order to use [CX](#), you must have openFrameworks installed. See <http://openframeworks.cc/download/> to download openFrameworks. Currently only version 0.8.0 of openFrameworks is supported by [CX](#).

Once you have installed openFrameworks, you can install [CX](#) by putting the contents of the [CX](#) repository into a sub-directory under openFrameworksDirectory/addons (typically openFrameworksDirectory/addons/ofxCX), where openFrameworksDirectory is where you put openFrameworks when you installed it.

To use [CX](#) in a project, use the openFrameworks project generator and select ofxCX as an addon. The project generator can be found in openFrameworksDirectory/projectGenerator.

In order to use the examples, do the following:

1. Use the oF project generator (in openFrameworksDirectory/projectGenerator) to create a new project that uses the ofxCX addon.

2. Go to the newly-created project directory (that you chose when creating the project in step 1) and go into the src subdirectory.
3. Delete all of the files in the src directory (main.cpp, testApp.h, and testApp.cpp). 4a. Copy the example .cpp file into this directory. 4b. If the example has a data folder, copy the contents of that folder into yourProject-Directory/bin/data. bin/data folders probably won't exist at this point. You can create them.
4. This step depends on your compiler, but you'll need to tell it to use the example source file that you copied in step 4a when it compiles the project (and possibly to specifically not use the files you deleted from the src directory in step 3).
5. Compile and run the project.

1.2 Examples and Tutorials

There are several examples that serve as tutorials for [CX](#). Some of the examples are on a specific topic and others are sample experiments that integrate together different features of [CX](#). The example files can be found in the [CX](#) directory (see [Installation](#)) in subfolders with names beginning with "example-".

Tutorials:

- soundObject - Tutorial covering a number of things that you can do with CX_SoundObjects, including loading sound files, combining sounds, and playing them.
- dataFrame - Tutorial covering use of CX_DataFrame, which is a container for storing data that is collected in an experiment.
- logging - Tutorial explaining how the error logging system of [CX](#) works and how you can use it in your experiments.

Experiments:

- basicChangeDetection - A very straightforward change-detection task demonstrating some of the features of [CX](#) like presentation of time-locked stimuli, keyboard response collection, and use of the CX_RandomNumber-Generator.
- advancedChangeDetection - This example expands on basicChangeDetection, including CX_DataFrame and CX_TrialController in order to simplify the experiment.
- nBack - Demonstrates advanced use of CX_SlidePresenter in the implementation of an N-Back task.

Misc.:

- helloWorld - A very basic getting started program.
- animation - A simple example of the most simple way to draw moving things in [CX](#). Also includes some mouse stuff: cursor movement, clicks, and scroll wheel activity.
- renderingTest - Includes several examples of how to draw stuff using ofPath (arbitrary lines), ofTexture (a kind of pixel buffer), ofImage (for opening image files: .png, .jpg, etc.), a variety of basic of drawing functions (ofCircle, ofRect, ofTriangle, etc.), and a number of [CX](#) drawing functions from the [CX::Draw](#) namespace.

2 Modules

3 Blocking Code

Blocking code is code that either takes a long time to complete or that waits until some event occurs before allowing code execution to continue. An example of blocking code that waits is

```
do {
    Input.pollEvents();
} while (Input.Keyboard.availableEvents() == 0);
```

This code waits until the keyboard has been used in some way. No code past it can be executed until the keyboard is used, which could take a long time. Any code that blocks while waiting for a human to do something is blocking.

An example of blocking code that takes a long time (or at least could take a long time) is

```
vector<double> d = CX::Util::sequence<double>(0, 1000000, .033);
```

which requires the allocation of about 300 MB of RAM. This code doesn't wait for anything to happen, it just takes a long time to execute.

Blocking code is bad because it prevents some parts of CX from working in some situations. It is not a cardinal sin and there are times when using blocking code is acceptable. However, blocking code should not be used when trying to present stimuli or when responses are being made. There is of course an exception to the responses rule, which is when your blocking code is explicitly polling for user input, e.g.:

```
while (!Input.pollEvents());
```

4 license

The code in this repository is available under the MIT License (https://secure.wikimedia.org/wikipedia/en/wiki/-Mit_license).

Copyright (c) 2014 Kyle Hardman

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

5 Module Index

5.1 Modules

Here is a list of all modules:

Data	8
Entry Point	9
Error Logging	10
Input Devices	11
Randomization	12
Sound	13
Timing	14
Utility	15
Video	16

6 Namespace Index

6.1 Namespace List

Here is a list of all documented namespaces with brief descriptions:

CX	17
CX::Draw	18
CX::Instances	20
CX::Private	20
CX::Util	21

7 Hierarchical Index

7.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

CX::CX_Clock	26
CX::Util::CX_CoordinateConverter	27
CX::CX_DataFrame	30
CX::CX_SafeDataFrame	67
CX::CX_DataFrameCell	37
CX::CX_DataFrameColumn	40
CX::CX_DataFrameRow	41

CX::Util::CX_DegreeToPixelConverter	41
CX::CX_Display	43
CX::CX_FinalSlideFunctionInfo_t	47
CX::CX_InputManager	48
CX::CX_Joystick	50
CX::CX_JoystickEvent_t	51
CX::CX_Keyboard	52
CX::CX_KeyEvent_t	53
CX::Util::CX_LengthToPixelConverter	54
CX::CX_Logger	55
CX::CX_Mouse	58
CX::CX_MouseEvent_t	59
CX::CX_RandomNumberGenerator	60
CX::CX_Slide_t	68
CX::CX_SlidePresenter	69
CX::CX_SlideTimingInfo_t	73
CX::CX_SoundObject	73
CX::CX_SoundObjectPlayer	77
CX::CX_SoundStream	79
CX::CX_SoundStreamConfiguration_t	83
CX::CX_SP_Configuration	84
CX::CX_SP_PresentationErrors_t	85
CX::Util::CX_TrialController	85

8 Class Index

8.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

CX::CX_Clock	26
CX::Util::CX_CoordinateConverter	27

CX::CX_DataFrame	30
CX::CX_DataFrameCell	37
CX::CX_DataFrameColumn	40
CX::CX_DataFrameRow	41
CX::Util::CX_DegreeToPixelConverter	41
CX::CX_Display	43
CX::CX_FinalSlideFunctionInfo_t	47
CX::CX_InputManager	48
CX::CX_Joystick	50
CX::CX_JoystickEvent_t	51
CX::CX_Keyboard	52
CX::CX_KeyEvent_t	53
CX::Util::CX_LengthToPixelConverter	54
CX::CX_Logger	55
CX::CX_Mouse	58
CX::CX_MouseEvent_t	59
CX::CX_RandomNumberGenerator	60
CX::CX_SafeDataFrame	67
CX::CX_Slide_t	68
CX::CX_SlidePresenter	69
CX::CX_SlideTimingInfo_t	73
CX::CX_SoundObject	73
CX::CX_SoundObjectPlayer	77
CX::CX_SoundStream	79
CX::CX_SoundStreamConfiguration_t	83
CX::CX_SP_Configuration	84
CX::CX_SP_PresentationErrors_t	85
CX::Util::CX_TrialController	85

9 File Index

9.1 File List

Here is a list of all documented files with brief descriptions:

advancedChangeDetectionTask.cpp	87
basicChangeDetectionTask.cpp	88
CX_Algorithm.h	??
CX_Clock.h	??
CX_ClockImplementations.h	??
CX_DataFrame.h	??
CX_DataFrame_attempts.h	??
CX_DataFrameCell.h	??
CX_Display.h	??
CX_Draw.h	??
CX_EntryPoint.h	??
CX_Events.h	??
CX_GLFWWindow_Compat.h	??
CX_InputManager.h	??
CX_Joystick.h	??
CX_Keyboard.h	??
CX_Logger.h	??
CX_LoggerChannel.h	??
CX_Mouse.h	??
CX_Private.h	??
CX_RandomNumberGenerator.h	??
CX_SlidePresenter.h	??
CX_SoundObject.h	??
CX_SoundObjectPlayer.h	??
CX_SoundStream.h	??
CX_SwappingThread.h	??

CX_TrialController.h	??
CX_TrialController_Class.h	??
CX_UnitConversion.h	??
CX_Uilities.h	??
myType.h	??
nBack.cpp	89

10 Module Documentation

10.1 Data

Classes

- class [CX::CX_DataFrame](#)
- class [CX::CX_DataFrameColumn](#)
- class [CX::CX_DataFrameRow](#)
- class [CX::CX_SafeDataFrame](#)
- class [CX::CX_DataFrameCell](#)

10.1.1 Detailed Description

This module is related to storing experimental data. [CX_DataFrame](#) is the most important class in this module.

10.2 Entry Point

Functions

- void `runExperiment` (void)

Variables

- `CX_Display CX::Instances::Display`
- `CX_InputManager CX::Instances::Input`
- `CX_Logger CX::Instances::Log`
- `CX_RandomNumberGenerator CX::Instances::RNG`

10.2.1 Detailed Description

The entry point provides access to a number of instances of classes that can be used by user code. It also provides declarations (but not definitions) of two functions which the user should define (`setupExperiment()` and `updateExperiment()`).

10.2.2 Function Documentation

10.2.2.1 `runExperiment (void)`

The user code should define a function with this name and type signature (takes no arguments and returns nothing). The user function will be called once setup is done for `CX`. When `runExperiment` returns, the program will exit.

```
void runExperiment (void) {  
    //Do your experiment.  
  
    return; //Return when done to exit the program. You don't have to explicitly return; you can just fall  
           off the end of the function.  
           //You can alternately call std::exit() or ofExit() at any point.  
}
```

10.2.3 Variable Documentation

10.2.3.1 `CX::CX_Display CX::Instances::Display`

An instance of `CX::CX_Display` that is lightly hooked into the `CX` backend. `setup()` is called for `Display` before `runExperiment()` is called.

10.2.3.2 `CX::CX_InputManager CX::Instances::Input`

An instance of `CX_InputManager` that is very lightly hooked into the `CX` backend.

10.2.3.3 `CX_Logger CX::Instances::Log`

This is an instance of `CX::CX_Logger` that is hooked into the `CX` backend. All log messages generated by `CX` and `openFrameworks` go through this instance.

10.2.3.4 `CX::CX_RandomNumberGenerator CX::Instances::RNG`

An instance of `CX_RandomNumberGenerator` that is (lightly) hooked into the `CX` backend.

10.3 Error Logging

Classes

- class [CX::CX_Logger](#)

Enumerations

- enum [CX::CX_LogLevel](#) {
 LOG_ALL, **LOG_VERBOSE**, **LOG_NOTICE**, **LOG_WARNING**,
 LOG_ERROR, **LOG_FATAL_ERROR**, **LOG_NONE** }

10.3.1 Detailed Description

10.3.2 Enumeration Type Documentation

10.3.2.1 enum [CX::CX_LogLevel](#) [strong]

Log levels for log messages. Depending on the log level chosen, the name of the level will be printed before the message. Depending on the settings set using `level()`, `levelForConsole()`, or `levelForFile()`, if the log level of a message is below the level set for the module or logging target it will not be printed. For example, if `LOG_ERROR` is the level for the console and `LOG_NOTICE` is the level for the module "test", then messages logged to the "test" module will be completely ignored if at verbose level (because of the module setting) and will not be printed to the console if they are below the level of an error (because of the console setting).

10.4 Input Devices

Classes

- class [CX::CX_InputManager](#)
- struct [CX::CX_JoystickEvent_t](#)
- class [CX::CX_Joystick](#)
- struct [CX::CX_KeyEvent_t](#)
- class [CX::CX_Keyboard](#)
- struct [CX::CX_MouseEvent_t](#)
- class [CX::CX_Mouse](#)

10.4.1 Detailed Description

There are a number of different classes that together perform the input handling functions of [CX](#). Start by looking at [CX::CX_InputManager](#) and the instance of that class that is created for you: [CX::Instances::Input](#).

For interfacing with serial ports, use `ofSerial` (<http://www.openframeworks.cc/documentation/communication/of-Serial.html>).

See Also

[CX::CX_InputManager](#) for the primary interface to input devices.

[CX::CX_Keyboard](#) for keyboard specific information.

[CX::CX_Mouse](#) for mouse specific information.

[CX::CX_Joystick](#) for joystick specific information.

10.5 Randomization

Classes

- class [CX::CX_RandomNumberGenerator](#)

10.5.1 Detailed Description

This module provides a class that is used for random number generation.

10.6 Sound

Classes

- struct [CX::CX_SoundStreamConfiguration_t](#)
- class [CX::CX_SoundObject](#)
- class [CX::CX_SoundObjectPlayer](#)
- class [CX::CX_SoundStream](#)

10.6.1 Detailed Description

10.7 Timing

Classes

- class [CX::CX_Clock](#)

Variables

- [CX_Clock CX::Instances::Clock](#)

10.7.1 Detailed Description

This module provides methods for timestamping events in experiments.

10.7.2 Variable Documentation

10.7.2.1 [CX_Clock CX::Instances::Clock](#)

An instance of [CX::CX_Clock](#) that is hooked into the [CX](#) backend. Anything in [CX](#) that requires timing information uses this instance. You should use this instance in your code and not make your own instance of [CX_Clock](#). You should never need another instance. You should never use another instance, as the experiment start times will not agree between instances.

10.8 Utility

Namespaces

- [CX::Util](#)

Classes

- class [CX::Util::CX_TrialController](#)
- class [CX::Util::CX_DegreeToPixelConverter](#)
- class [CX::Util::CX_LengthToPixelConverter](#)
- class [CX::Util::CX_CoordinateConverter](#)

10.8.1 Detailed Description

10.9 Video

Namespaces

- [CX::Draw](#)

Classes

- class [CX::CX_Display](#)
- struct [CX::CX_FinalSlideFunctionInfo_t](#)
- struct [CX::CX_SP_PresentationErrors_t](#)
- struct [CX::CX_SP_Configuration](#)
- struct [CX::CX_SlideTimingInfo_t](#)
- struct [CX::CX_Slide_t](#)
- class [CX::CX_SlidePresenter](#)

10.9.1 Detailed Description

This module is related to creating and presenting visual stimuli. Mostly, it is responsible for controlling presentation timing of stimuli rather than actually creating the stimuli.

Almost all of the actual drawing of stimuli is done using openFrameworks functions. A lot of the common functions can be found in ofGraphics.h (<http://www.openframeworks.cc/documentation/graphics/of-Graphics.html>), but there are a lot of other ways to draw stimuli: see the graphics and 3d sections of this page: <http://www.openframeworks.cc/documentation/>.

11 Namespace Documentation

11.1 CX Namespace Reference

Namespaces

- [Draw](#)
- [Instances](#)
- [Private](#)
- [Util](#)

Classes

- class [CX_Clock](#)
- class [CX_DataFrame](#)
- class [CX_DataFrameColumn](#)
- class [CX_DataFrameRow](#)
- class [CX_SafeDataFrame](#)
- class [CX_DataFrameCell](#)
- class [CX_Display](#)
- class [CX_InputManager](#)
- struct [CX_JoystickEvent_t](#)
- class [CX_Joystick](#)
- struct [CX_KeyEvent_t](#)
- class [CX_Keyboard](#)
- class [CX_Logger](#)
- struct [CX_MouseEvent_t](#)
- class [CX_Mouse](#)
- class [CX_RandomNumberGenerator](#)
- struct [CX_FinalSlideFunctionInfo_t](#)
- struct [CX_SP_PresentationErrors_t](#)
- struct [CX_SP_Configuration](#)
- struct [CX_SlideTimingInfo_t](#)
- struct [CX_Slide_t](#)
- class [CX_SlidePresenter](#)
- class [CX_SoundObject](#)
- class [CX_SoundObjectPlayer](#)
- struct [CX_SoundStreamConfiguration_t](#)
- class [CX_SoundStream](#)

Typedefs

- typedef long long [CX_Micros](#)
- typedef int64_t [CX_RandomInt_t](#)
- typedef [CX_SoundStreamConfiguration_t](#) [CX_SoundObjectPlayerConfiguration_t](#)

Enumerations

- enum [CX_LogLevel](#) {
 LOG_ALL, **LOG_VERBOSE**, **LOG_NOTICE**, **LOG_WARNING**,
 LOG_ERROR, **LOG_FATAL_ERROR**, **LOG_NONE** }
- enum **CX_MouseButtons** : int { **LEFT_MOUSE** = OF_MOUSE_BUTTON_LEFT, **MIDDLE_MOUSE** = OF_MOUSE_BUTTON_MIDDLE, **RIGHT_MOUSE** = OF_MOUSE_BUTTON_RIGHT }
- enum [CX_SP_ErrorMode](#) { [CX_SP_ErrorMode::PROPAGATE_DELAYS](#), [CX_SP_ErrorMode::FIX_TIMING_FROM_FIRST_SLIDE](#) }

Functions

- std::ostream & **operator**<< (std::ostream &os, const [CX_DataFrameCell](#) &cell)
- std::ostream & **operator**<< (std::ostream &os, const [CX_JoystickEvent_t](#) &ev)
- std::istream & **operator**>> (std::istream &is, [CX_JoystickEvent_t](#) &ev)
- std::ostream & **operator**<< (std::ostream &os, const [CX_KeyEvent_t](#) &ev)
- std::istream & **operator**>> (std::istream &is, [CX_KeyEvent_t](#) &ev)
- std::ostream & **operator**<< (std::ostream &os, const [CX_MouseEvent_t](#) &ev)
- std::istream & **operator**>> (std::istream &is, [CX_MouseEvent_t](#) &ev)

11.1.1 Detailed Description

This namespace contains all of the symbols related to [CX](#).

11.1.2 Enumeration Type Documentation

11.1.2.1 enum **CX::CX_SP_ErrorMode** [strong]

The settings in this enum are related to what a [CX_SlidePresenter](#) does when it encounters a timing error. Timing errors are probably almost exclusively related to one slide being presented for too long.

The PROPAGATE_DELAYS setting causes the slide presenter to handle these errors by moving the start time of all future stimuli back by the number of extra frame that the erroneous slide used. This makes the durations of all future stimuli correct, so that there is only an error in the duration of one slide.

An alternative option is to try to keep the onsets of all slides as constant as possible relative to each other. This means that if one slide is presented for an extra frame, the next slide will be presented for one frame less than it should have been. If one slide is presented for several extra frames (this should almost never happen), the next slide may be skipped altogether. However, this mode ([FIX_TIMING_FROM_FIRST_SLIDE](#)) does not completely work currently so it should not be used.

Enumerator

PROPAGATE_DELAYS This mode handles timing errors by changing the onset times of future stimuli so.

FIX_TIMING_FROM_FIRST_SLIDE This does not work currently.

11.2 CX::Draw Namespace Reference

Functions

- ofPath [squircleToPath](#) (double radius, double amount=0.9)
- ofPath [starToPath](#) (int numberOfPoints, double innerRadius, double outerRadius)

- void [star](#) (ofPoint center, int numberOfPoints, float innerRadius, float outerRadius, ofColor lineColor, ofColor fillColor, float lineWidth=1, float rotationDeg=0)
- void [centeredString](#) (int x, int y, std::string s, ofTrueTypeFont &font)
- void **centeredString** (ofPoint center, std::string s, ofTrueTypeFont &font)
- ofPixels **greyscalePattern** (const CX_PatternProperties_t &patternProperties)
- ofPixels **gaborToPixels** (const CX_GaborProperties_t &properties)
- ofTexture **gaborToTexture** (const CX_GaborProperties_t &properties)
- void **gabor** (int x, int y, const CX_GaborProperties_t &properties)

11.2.1 Detailed Description

This namespace contains functions for drawing certain complex stimuli.

This namespace contains a variety of [CX](#) drawing functions.

11.2.2 Function Documentation

11.2.2.1 void CX::Draw::centeredString (int x, int y, std::string s, ofTrueTypeFont & font)

Draws a string centered on a given location using the given font. Strings are normally drawn such that the x coordinate gives the left edge of the string and the y coordinate gives a line on which the letters will be drawn, where some characters (like y or g) can descend below the line.

Parameters

<i>x</i>	The x coordinate of the center of the string.
<i>y</i>	The y coordinate of the center of the string.
<i>s</i>	The string to draw.
<i>font</i>	A font that has already been prepared for use.

11.2.2.2 ofPath CX::Draw::squircleToPath (double radius, double amount = 0.9)

This function draws an approximation of a squircle (<http://en.wikipedia.org/wiki/Squircle>) using Bezier curves. The squircle will be centered on (0,0) in the ofPath.

Parameters

<i>radius</i>	The radius of the largest circle that can be enclosed in the squircle.
<i>amount</i>	The "squirliness" of the squircle. The default (0.9) seems like a pretty good amount for a good approximation of a squircle, but different amounts can give different sorts of shapes.

Returns

An ofPath containing the squircle.

11.2.2.3 void CX::Draw::star (ofPoint center, int numberOfPoints, float innerRadius, float outerRadius, ofColor lineColor, ofColor fillColor, float lineWidth = 1, float rotationDeg = 0)

This draws an N-pointed star.

Parameters

<i>center</i>	The point at the center of the star.
<i>numberOfPoints</i>	The number of points in the star.
<i>innerRadius</i>	The distance from the center of the star to where the inner points of the star hit.
<i>outerRadius</i>	The distance from the center of the star to the outer points of the star.
<i>lineColor</i>	The color of the lines going around the edge of the star.
<i>fillColor</i>	The color used to fill in the center of the star.
<i>lineWidth</i>	The width of the lines.
<i>rotationDeg</i>	The number of degrees to rotate the star. 0 degrees has one point of the star pointing up. Positive values rotate the star counter-clockwise.

Returns

An ofPath containing the star.

11.2.2.4 ofPath CX::Draw::starToPath (int *numberOfPoints*, double *innerRadius*, double *outerRadius*)

This draws an N-pointed star to an ofPath. The star will be centered on (0,0) in the ofPath.

Parameters

<i>numberOfPoints</i>	The number of points in the star.
<i>innerRadius</i>	The distance from the center of the star at which the inner points of the star hit.
<i>outerRadius</i>	The distance from the center of the star to the outer points of the star.

Returns

An ofPath containing the star.

11.3 CX::Instances Namespace Reference

Variables

- [CX_Clock](#) Clock
- [CX_Display](#) Display
- [CX_InputManager](#) Input
- [CX_Logger](#) Log
- [CX_RandomNumberGenerator](#) RNG

11.3.1 Detailed Description

This namespace contains instances of some classes that are fundamental to the functioning of [CX](#).

11.4 CX::Private Namespace Reference

Functions

- void **glfwErrorCallback** (int code, const char *message)
- CX_Events & **getEvents** (void)
- CX_GLVersion **getOpenGLVersion** (void)
- CX_GLVersion **getGLSLVersion** (void)

- bool **glFenceSyncSupported** (void)
- bool **glVersionAtLeast** (int desiredMajor, int desiredMinor, int desiredRelease=0)
- int **glCompareVersions** (CX_GLVersion a, CX_GLVersion b)
- CX_GLVersion **getGLSLVersionFromGLVersion** (CX_GLVersion glVersion)

Variables

- GLFWwindow * **glfwContext** = NULL
- ofPtr
 < ofAppGLFWCompatibilityWindow > **window**

11.4.1 Detailed Description

This namespace contains symbols that may be visible in user code but which should not be used by user code.

11.5 CX::Util Namespace Reference

Classes

- class [CX_TrialController](#)
- class [CX_DegreeToPixelConverter](#)
- class [CX_LengthToPixelConverter](#)
- class [CX_CoordinateConverter](#)

Enumerations

- enum **CX_RoundingConfiguration** { **ROUND_TO_NEAREST**, **ROUND_UP**, **ROUND_DOWN**, **ROUND_TOWARD_ZERO** }

Functions

- float **degreesToPixels** (float degrees, float pixelsPerUnit, float viewingDistance)
- float [pixelsToDegrees](#) (float pixels, float pixelsPerUnit, float viewingDistance)
- bool [checkOFVersion](#) (int versionMajor, int versionMinor, int versionPatch)
- int **getSampleCount** (void)
- template<typename T >
 std::vector< T > [arrayToVector](#) (T arr[], unsigned int arraySize)
- template<typename T >
 std::vector< T > [sequence](#) (T start, T end, T stepSize)
- template<typename T >
 std::vector< T > [sequenceSteps](#) (T start, unsigned int steps, T stepSize)
- template<typename T >
 std::vector< T > [sequenceAlong](#) (T start, T end, unsigned int steps)
- template<typename T >
 std::vector< T > [intVector](#) (T start, T end)
- template<typename T >
 std::vector< T > [repeat](#) (T value, unsigned int times)
- template<typename T >
 std::vector< T > [repeat](#) (std::vector< T > values, unsigned int times, unsigned int each=1)

- `template<typename T >`
`std::vector< T > repeat (std::vector< T > values, std::vector< unsigned int > each, unsigned int times=1)`
- `template<typename T >`
`std::string vectorToString (std::vector< T > values, std::string delimiter=",", int significantDigits=8)`
- `bool writeToFile (std::string filename, std::string data, bool append=true)`
- `double round (double d, int roundingPower, CX_RoundingConfiguration c)`

11.5.1 Detailed Description

This namespace contains a variety of utility functions.

11.5.2 Function Documentation

11.5.2.1 `template<typename T > std::vector< T > CX::Util::arrayToVector (T arr[], unsigned int arraySize)`

Copies `arraySize` elements of an array of `T` to a `vector<T>`.

Template Parameters

<code>< T ></code>	The type of the array. Is often inferred by the compiler.
--------------------------	---

Parameters

<code>arr</code>	The array of data to put into the vector.
<code>arraySize</code>	The length of the array, or the number of elements to copy from the array if not all of the elements are wanted.

Returns

The elements in a vector.

11.5.2.2 `bool CX::Util::checkOFVersion (int versionMajor, int versionMinor, int versionPatch)`

Checks that the version of oF that is used during compilation matches the requested version. If the desired version was 0.7.1, simply input (0, 7, 1) as the arguments. A warning will be logged if the versions don't match.

Returns

True if the versions match, false otherwise.

11.5.2.3 `template<typename T > std::vector< T > CX::Util::intVector (T start, T end)`

Creates a vector of integers going from `start` to `end`. `start` may be greater than `end`, in which case the returned values will be in descending order. This is similar to using `CX::sequence`, but the step size is fixed to 1 and it works properly when trying to create a descending sequence of unsigned integers.

Returns

A vector of the values in the sequence.

11.5.2.4 `float CX::Util::pixelsToDegrees (float pixels, float pixelsPerUnit, float viewingDistance)`

The inverse of `CX::Util::degreesToPixels()`.

11.5.2.5 `template<typename T> std::vector< T > CX::Util::repeat (T value, unsigned int times)`

Repeats value "times" times.

Parameters

<i>value</i>	The value to be repeated.
<i>times</i>	The number of times to repeat the value.

Returns

A vector containing times copies of the repeated value.

11.5.2.6 `template<typename T> std::vector< T > CX::Util::repeat (std::vector< T > values, unsigned int times, unsigned int each = 1)`

Repeats the elements of values. Each element of values is repeated "each" times and then the process of repeating the elements is repeated "times" times.

Parameters

<i>values</i>	Vector of values to be repeated.
<i>times</i>	The number of times the process should be performed.
<i>each</i>	Number of times each element of values should be repeated.

Returns

A vector of the repeated values.

11.5.2.7 `template<typename T> std::vector< T > CX::Util::repeat (std::vector< T > values, std::vector< unsigned int > each, unsigned int times = 1)`

Repeats the elements of values. Each element of values is repeated "each" times and then the process of repeating the elements is repeated "times" times.

Parameters

<i>values</i>	Vector of values to be repeated.
<i>each</i>	Number of times each element of values should be repeated. Must be the same length as values. If not, an error is logged and an empty vector is returned.
<i>times</i>	The number of times the process should be performed.

Returns

A vector of the repeated values.

11.5.2.8 `double CX::Util::round (double d, int roundingPower, CX_RoundingConfiguration c)`

Rounds the given double to the given power of 10.

Parameters

<i>d</i>	The number to be rounded.
<i>roundingPower</i>	The power of 10 to round d to. For example, if roundingPower is 0, d is rounded to the one's place ($10^0 == 1$). If roundingPower is -3, d is rounded to the thousandth's place ($10^{-3} = .001$). If roundingPower is 1, d is rounded to the ten's place.

<i>c</i>	The type of rounding to do, from the CX::Util::CX_RoundingConfiguration enum. You can round up, down, to nearest, and toward zero.
----------	--

Returns

The rounded value.

11.5.2.9 `template<typename T> std::vector< T > CX::Util::sequence (T start, T end, T stepSize)`

Creates a sequence of numbers from start to end by steps of size stepSize. start may be geater than end, but only if stepSize is less than 0. If start is less than end, stepSize must be greater than 0.

Example call: `sequence<double>(1, 3.3, 2)` results in {1, 3}

Parameters

<i>start</i>	The start of the sequence. You are guaranteed to get this value in the sequence.
<i>end</i>	The number past which the sequence should end. You are not guaranteed to get this value.
<i>stepSize</i>	A nonzero number.

11.5.2.10 `template<typename T> std::vector< T > CX::Util::sequenceAlong (T start, T end, unsigned int outputLength)`

Creates a sequence from start to end, where the size of each step is chosen so that the length of the sequence if equal to outputLength.

Parameters

<i>start</i>	The value at which to start the sequence.
<i>end</i>	The value to which to end the sequence.
<i>outputLength</i>	The number of elements in the returned sequence.

Returns

A vector containing the sequence.

11.5.2.11 `template<typename T> std::vector< T > CX::Util::sequenceSteps (T start, unsigned int steps, T stepSize)`

Make a sequence starting from start and taking steps steps of stepSize.

`sequenceSteps(1.5, 4, 2.5);`

Creates the sequence {1.5, 4, 6.5, 9, 11.5}

Parameters

<i>start</i>	Value from which to start.
<i>steps</i>	The number of steps to take.
<i>stepSize</i>	The size of each step.

Returns

A vector containing the sequence.

11.5.2.12 `bool CX::Util::writeToFile (std::string filename, std::string data, bool append = true)`

Writes data to a file, either appending the data to an existing file or creating a new file, overwriting any existing file with the given filename.

Parameters

<i>filename</i>	Name of the file to write to. If it is a relative file name, it will be placed relative the the data directory.
<i>data</i>	The data to write
<i>append</i>	If true, data will be appended to an existing file, if it exists. If append is false, any existing file will be overwritten and a warning will be logged. If no file exists, a new one will be created.

Returns

True if an error was encountered while writing the file, true otherwise. If there was an error, an error message will be logged.

12 Class Documentation

12.1 CX::CX_Clock Class Reference

```
#include <CX_Clock.h>
```

Public Types

- typedef std::chrono::steady_clock **CX_InternalClockType**

Public Member Functions

- void **precisionTest** (unsigned int iterations)
- CX_Micros [getTime](#) (void)
- CX_Micros [getSystemTime](#) (void)
- CX_Micros [getExperimentStartTime](#) (void)
- std::string [getExperimentStartDateTimeString](#) (std::string format="%Y-%b-%e %h-%M-%S %a")

Static Public Member Functions

- static std::string [getDateTimeString](#) (std::string format="%Y-%b-%e %h-%M-%S %a")

12.1.1 Detailed Description

This class is responsible for getting timestamps for anything requiring timestamps. The way to get timing information is the function [getTime\(\)](#). It returns the current time relative to the start of the experiment in microseconds (on most systems, see [getTickPeriod\(\)](#) to check the actual precision).

An instance of this class is preinstantiated for you. See [CX::Instances::Clock](#).

12.1.2 Member Function Documentation

12.1.2.1 std::string CX_Clock::getDateTimeString (std::string *format* = "%Y-%b-%e %h-%M-%S %a") [static]

This function returns a string containing the local time encoded according to some format.

Parameters

<i>format</i>	See http://pocoproject.org/docs/Poco.DateTimeFormatter.html#4684 for documentation of the format. E.g. "%Y/%m/%d %H:%M:%S" gives "year/month/day 24Hour-Clock:minute:second" with some zero-padding for most things. The default "%Y-%b-%e %h-%M-%S %a" is "yearWithCentury-abbreviatedMonthName-nonZeroPaddedDay 12HourClock-minuteZeroPadded-secondZeroPadded am/pm".
---------------	--

12.1.2.2 std::string CX_Clock::getExperimentStartDateTimeString (std::string *format* = "%Y-%b-%e %h-%M-%S %a")

Get a string representing the date/time of the start of the experiment encoded according to a format.

Parameters

<i>format</i>	See getDateTimeString() for the definition of the format.
---------------	---

12.1.2.3 CX_Micros CX_Clock::getExperimentStartTime (void)

Get the start time of the experiment in system time. The returned value can be compared with the result of [getSystemTime\(\)](#).

12.1.2.4 CX_Micros CX_Clock::getSystemTime (void)

This function returns the current system time in microseconds.

This cannot be converted to time/day in any meaningful way. Use [getDateTimeString\(\)](#) for that.

Returns

A time value that can be compared to the result of other calls to this function and to [getExperimentStartTime\(\)](#).

12.1.2.5 CX_Micros CX_Clock::getTime (void)

This function returns the current time relative to the start of the experiment in microseconds. The start of the experiment is defined by default as when the [CX_Clock](#) instance named Clock (instantiated in this file) is constructed (typically the beginning of program execution).

The documentation for this class was generated from the following files:

- CX_Clock.h
- CX_Clock.cpp

12.2 CX::Util::CX_CoordinateConverter Class Reference

```
#include <CX_UnitConversion.h>
```

Public Member Functions

- [CX_CoordinateConverter](#) (ofPoint origin, bool invertX, bool invertY, bool invertZ=false)
- ofPoint [operator\(\)](#) (ofPoint p)
- ofPoint [operator\(\)](#) (float x, float y, float z=0)
- void [setAxisInversion](#) (bool invertX, bool invertY, bool invertZ=false)
- void [setOrigin](#) (ofPoint newOrigin)
- void [setUnitConverter](#) (CX_BaseUnitConverter *converter)

12.2.1 Detailed Description

This helper class is used for converting from a somewhat user-defined coordinate system into the standard computer monitor coordinate system. When user coordinates are input into this class, they will be converted into the standard monitor coordinate system. This lets you use this class to allow you to use coordinates in your own system and convert those coordinates into the standard coordinates that are used by the drawing functions of openFrameworks.

See [setUnitConverter\(\)](#) for a way to do change the units of the coordinate system to, for example, inches or degrees of visual angle.

Example use:

```
CX_CoordinateConverter conv(Display.  
    getCenterOfDisplay(), false, true); //Make the center of the display the origin and  
    invert  
//the Y-axis. This makes positive x values go to the right and positive y values go up from the center of  
    the display.  
ofSetColor(255, 0, 0); //Draw a red circle in the center of the display.  
ofCircle(conv(0, 0), 20);  
ofSetColor(0, 255, 0); //Draw a green circle 100 pixels to the right of the center.  
ofCircle(conv(100, 0), 20);  
ofSetColor(0, 0, 255); //Draw a blue circle 100 pixels above the center (inverted y-axis).  
ofCircle(conv(0, 100), 20);
```

12.2.2 Constructor & Destructor Documentation

12.2.2.1 CX::Util::CX_CoordinateConverter::CX_CoordinateConverter (ofPoint *origin*, bool *invertX*, bool *invertY*, bool *invertZ* = false)

Constructs a coordinate converter with the given settings.

Parameters

<i>origin</i>	The location within the standard coordinate system at which the origin (the point at which the x, y, and z values are 0) of the user-defined coordinate system is located. If, for example, you want the center of the display to be the origin within your user-defined coordinate system, you could use CX_Display::getCenterOfDisplay() as the value for this argument.
<i>invertX</i>	Invert the x-axis from the default, which is that x increases to the right.
<i>invertY</i>	Invert the y-axis from the default, which is that y increases downward.
<i>invertZ</i>	Invert the z-axis from the default, which is that z increases toward the user (i.e. pointing out of the front of the screen).

12.2.3 Member Function Documentation

12.2.3.1 ofPoint CX::Util::CX_CoordinateConverter::operator() (ofPoint *p*)

The primary method of conversion between coordinate systems.

Parameters

<i>p</i>	The point in user coordinates that should be converted to standard coordinates.
----------	---

Returns

The point in standard coordinates.

Example use:

```
CX_CoordinateConverter cc(ofPoint(200,200), false, true);
```

```
ofPoint p(-50, 100); //P is in user-defined coordinates, 50 units left and 100 units above the origin.
ofPoint res = cc(p); //Use operator() to convert from the user system to the standard system.
//res should contain (150, 100) due to the inverted y axis.
```

12.2.3.2 ofPoint CX::Util::CX_CoordinateConverter::operator() (float x, float y, float z = 0)

Equivalent to a call to operator()(ofPoint(x, y, z)).

12.2.3.3 void CX::Util::CX_CoordinateConverter::setAxisInversion (bool invertX, bool invertY, bool invertZ = false)

Sets whether each axis within the user-defined system is inverted from the standard coordinate system.

Parameters

<i>invertX</i>	Invert the x-axis from the default, which is that x increases to the right.
<i>invertY</i>	Invert the y-axis from the default, which is that y increases downward.
<i>invertZ</i>	Invert the z-axis from the default, which is that z increases toward the user (i.e. pointing out of the front of the screen).

12.2.3.4 void CX::Util::CX_CoordinateConverter::setOrigin (ofPoint newOrigin)

Sets the location within the standard coordinate system at which the origin of the user-defined coordinate system is located.

Parameters

<i>newOrigin</i>	The location within the standard coordinate system at which the origin (the point at which the x, y, and z values are 0) of the user-defined coordinate system is located. If, for example, you want the center of the display to be the origin within your user-defined coordinate system, you could use CX_Display::getCenterOfDisplay() as the value for this argument.
------------------	--

12.2.3.5 void CX::Util::CX_CoordinateConverter::setUnitConverter (CX_BaseUnitConverter * converter)

Sets the unit converter that will be used when converting the coordinate system. In this way you can convert both the coordinate system in use and the units used by the coordinate system. See [CX_DegreeToPixelConverter](#) and [CX_LengthToPixelConverter](#) for examples of the converters that can be used.

Example use:

```
//At global scope:
CX_CoordinateConverter conv(Display.
    getCenterOfDisplay(), false, true);
CX_DegreeToPixelConverter d2p(35, 70);

//During setup:
conv.setUnitConverter(&d2p); //Use degrees of visual angle as the units of the user coordinate system.

//Draw a blue circle 2 degrees of visual angle to the left of the origin and 3 degrees above the origin
    (inverted y-axis).
ofSetColor(0, 0, 255);
ofCircle(conv(-2, 3), 20);
```

Parameters

<i>converter</i>	A pointer to an instance of a class that is a CX_BaseUnitConverter or which has inherited from that class. See CX_UnitConversion.h/cpp for the implementation of CX_LengthToPixelConverter to see an example of how to create you own converter.
------------------	--

Note

The origin of the coordinate converter must be in the units that result from the unit conversion. E.g. if you are converting the units from degrees to pixels, the origin must be in pixels. See [setOrigin\(\)](#).

The unit converter passed to this function must continue to exist throughout the lifetime of the coordinate converter. It is not copied.

The documentation for this class was generated from the following files:

- CX_UnitConversion.h
- CX_UnitConversion.cpp

12.3 CX::CX_DataFrame Class Reference

```
#include <CX_DataFrame.h>
```

Inherited by [CX::CX_SafeDataFrame](#) [protected].

Public Types

- typedef std::vector
 < [CX_DataFrameCell](#) >
 ::size_type **rowIndex_t**

Public Member Functions

- [CX_DataFrame](#) & **operator=** ([CX_DataFrame](#) &df)
- [CX_DataFrameCell](#) **operator()** (std::string column, [rowIndex_t](#) row)
- [CX_DataFrameCell](#) **operator()** ([rowIndex_t](#) row, std::string column)
- [CX_DataFrameCell](#) **at** ([rowIndex_t](#) row, std::string column)
- [CX_DataFrameCell](#) **at** (std::string column, [rowIndex_t](#) row)
- [CX_DataFrameColumn](#) **operator[]** (std::string column)
- [CX_DataFrameRow](#) **operator[]** ([rowIndex_t](#) row)
- void **appendRow** ([CX_DataFrameRow](#) row)
- std::string **print** (std::string delimiter="\t", bool printRowNumbers=true)
- std::string **print** (const std::set< std::string > &columns, std::string delimiter="\t", bool printRowNumbers=true)
- std::string **print** (const std::vector< [rowIndex_t](#) > &rows, std::string delimiter="\t", bool printRowNumbers=true)
- std::string **print** (const std::set< std::string > &columns, const std::vector< [rowIndex_t](#) > &rows, std::string delimiter="\t", bool printRowNumbers=true)
- bool **printToFile** (std::string filename, std::string delimiter="\t", bool printRowNumbers=true)
- bool **printToFile** (std::string filename, const std::set< std::string > &columns, std::string delimiter="\t", bool printRowNumbers=true)
- bool **printToFile** (std::string filename, const std::vector< [rowIndex_t](#) > &rows, std::string delimiter="\t", bool printRowNumbers=true)
- bool **printToFile** (std::string filename, const std::set< std::string > &columns, const std::vector< [rowIndex_t](#) > &rows, std::string delimiter="\t", bool printRowNumbers=true)
- bool **readFromFile** (std::string filename, std::string cellDelimiter="\t", std::string vectorEncloser="")
- void **clear** (void)
- bool **deleteColumn** (std::string columnName)
- bool **deleteRow** ([rowIndex_t](#) row)
- std::vector< std::string > **columnNames** (void)

- `rowIndex_t` [getRowCount](#) (void)
Returns the number of rows in the data frame.
- `bool` [reorderRows](#) (const vector< CX_DataFrame::rowIndex_t > &newOrder)
- [CX_DataFrame copyRows](#) (vector< CX_DataFrame::rowIndex_t > rowOrder)
- [CX_DataFrame copyColumns](#) (vector< std::string > columns)
- `void` [shuffleRows](#) (void)
- `void` [shuffleRows](#) (CX_RandomNumberGenerator &rng)
- `template<typename T >`
`std::vector< T >` [copyColumn](#) (std::string column)
- `void` [setRowCount](#) (rowIndex_t rowCount)
- `void` [addColumn](#) (std::string columnName)

Protected Member Functions

- `void` [_resizeToFit](#) (std::string column, rowIndex_t row)
- `void` [_resizeToFit](#) (rowIndex_t row)
- `void` [_resizeToFit](#) (std::string column)
- `void` [_equalizeRowLengths](#) (void)

Protected Attributes

- `std::map< std::string, vector< CX_DataFrameCell > >` [_data](#)
- `rowIndex_t` [_rowCount](#)

Friends

- `class` [CX_DataFrameRow](#)
- `class` [CX_DataFrameColumn](#)

12.3.1 Detailed Description

This class provides an easy way to store data from an experiment and output that data to a file at the end of the experiment. A [CX_DataFrame](#) is a square two-dimensional array of cells, but each cell is capable of holding a vector of data. Each cell is indexed with a column name (a string) and a row number. Cells can store many different kinds of data and the data can be inserted or extracted easily. The standard method of storing data is to use `operator()`, which dynamically resizes the data frame. When an experimental session is complete, the data can be written to a file using [printToFile\(\)](#).

See the example `dataFrame.cpp` for thorough examples of how to use a [CX_DataFrame](#).

12.3.2 Member Function Documentation

12.3.2.1 `void CX_DataFrame::addColumn (std::string columnName)`

Adds a column to the data frame.

Parameters

<i>columnName</i>	The name of the column to add. If a column with that name already exists in the data frame, a warning will be logged.
-------------------	---

12.3.2.2 void CX_DataFrame::appendRow (CX_DataFrameRow row)

Appends the row to the end of the data frame.

Parameters

<i>row</i>	The row to add.
------------	-----------------

Note

If the row has columns that do not exist in the data frame, those columns will be added to the data frame.

12.3.2.3 CX_DataFrameCell CX_DataFrame::at (rowIndex_t row, std::string column)

Access the cell at the given row and column with bounds checking. Throws a std::exception and logs an error if either the row or column is out of bounds.

Parameters

<i>row</i>	The row number.
<i>column</i>	The column name.

Returns

A [CX_DataFrameCell](#) that can be read from or written to.

12.3.2.4 void CX_DataFrame::clear (void)

Deletes the contents of the data frame. Resizes the data frame to have no rows and no columns.

12.3.2.5 std::vector< std::string > CX_DataFrame::columnNames (void)

Returns a vector containing the names of the columns in the data frame.

Returns

Vector of strings with the column names.

12.3.2.6 template<typename T > std::vector< T > CX::CX_DataFrame::copyColumn (std::string column)

Makes a copy of the data contained in the named column, converting it to the specified type (such a conversion must be possible).

Parameters

<i>column</i>	The name of the column to copy data from.
---------------	---

Returns

A vector containing the copied data.

12.3.2.7 CX_DataFrame CX_DataFrame::copyColumns (vector< std::string > *columns*)

Copies the specified columns into a new data frame.

Parameters

<i>columns</i>	A vector of column names to copy out. If a requested column is not found, a warning will be logged, but the function will otherwise complete successfully.
----------------	--

Returns

A [CX_DataFrame](#) containing the specified columns.

12.3.2.8 CX_DataFrame CX_DataFrame::copyRows (vector< CX_DataFrame::rowIndex_t > rowOrder)

Creates [CX_DataFrame](#) containing a copy of the rows specified in rowOrder. The new data frame is not linked to the existing data frame.

Parameters

<i>rowOrder</i>	A vector of CX_DataFrame::rowIndex_t containing the rows from this data frame to be copied out. The indices in rowOrder may be in any order: They don't need to be ascending. Additionally, the same row to be copied may be specified multiple times.
-----------------	--

Returns

A [CX_DataFrame](#) containing the rows specified in rowOrder.

12.3.2.9 bool CX_DataFrame::deleteColumn (std::string columnName)

Deletes the given column of the data frame.

Parameters

<i>columnName</i>	The name of the column to delete. If the column is not in the data frame, a warning will be logged.
-------------------	---

Returns

True if the column was found and deleted, false if it was not found.

12.3.2.10 bool CX_DataFrame::deleteRow (rowIndex_t row)

Deletes the given row of the data frame.

Parameters

<i>row</i>	The row to delete (0 indexed). If row is greater than or equal to the number of rows in the data frame, a warning will be logged.
------------	---

Returns

True if the row was in bounds and was deleted, false if the row was out of bounds.

12.3.2.11 CX_DataFrameCell CX_DataFrame::operator() (std::string column, rowIndex_t row)

Access the cell at the given row and column. If the row or column is out of bounds, the data frame will be dynamically resized in order to fit the row or column.

Parameters

<i>row</i>	The row number.
<i>column</i>	The column name.

Returns

A [CX_DataFrameCell](#) that can be read from or written to.

12.3.2.12 CX_DataFrame & CX_DataFrame::operator= (CX_DataFrame & df)

Copy the contents of another [CX_DataFrame](#).

Parameters

<i>df</i>	The data frame to copy.
-----------	-------------------------

Returns

A reference to this data frame.

Note

The contents of this data frame are deleted during the copy.

12.3.2.13 std::string CX_DataFrame::print (std::string delimiter = "\t", bool printRowNumbers = true)

Reduced argument version of [print\(\)](#). Prints all rows and columns.

12.3.2.14 std::string CX_DataFrame::print (const std::set< std::string > & columns, std::string delimiter = "\t", bool printRowNumbers = true)

Reduced argument version of [print\(\)](#). Prints all rows and the selected columns.

12.3.2.15 std::string CX_DataFrame::print (const std::vector< rowIndex_t > & rows, std::string delimiter = "\t", bool printRowNumbers = true)

Reduced argument version of [print\(\)](#). Prints all columns and the selected rows.

12.3.2.16 std::string CX_DataFrame::print (const std::set< std::string > & columns, const std::vector< rowIndex_t > & rows, std::string delimiter = "\t", bool printRowNumbers = true)

Prints the selected rows and columns of the data frame to a string. Each cell of the data frame will be separated with the selected delimiter.

Each row of the data frame will be ended with a new line (whatever `std::endl` evaluates to, typically `"\r\n"`).

Parameters

<i>columns</i>	Columns to print. Column names not found in the data frame will be ignored with a warning.
<i>rows</i>	Rows to print. Row indices not found in the data frame will be ignored with a warning.
<i>delimiter</i>	Delimiter to be used between cells of the data frame. Using comma or semicolon for the delimiter is not recommended because semicolons are used as element delimiters in the string-encoded vectors stored in the data frame and commas are used for element delimiters within each element of the string-encoded vectors.
<i>printRowNumbers</i>	If true, a column will be printed with the header "rowNumber" with the contents of the column being the selected row indices. If false, no row numbers will be printed.

Returns

A string containing the printed version of the data frame.

12.3.2.17 `bool CX_DataFrame::printToFile (std::string filename, std::string delimiter = "\t", bool printRowNumbers = true)`

Reduced argument version of [printToFile\(\)](#). Prints all rows and columns.

12.3.2.18 `bool CX_DataFrame::printToFile (std::string filename, const std::set< std::string > & columns, std::string delimiter = "\t", bool printRowNumbers = true)`

Reduced argument version of [printToFile\(\)](#). Prints all rows and the selected columns.

12.3.2.19 `bool CX_DataFrame::printToFile (std::string filename, const std::vector< rowIndex_t > & rows, std::string delimiter = "\t", bool printRowNumbers = true)`

Reduced argument version of [printToFile\(\)](#). Prints all columns and the selected rows.

12.3.2.20 `bool CX_DataFrame::printToFile (std::string filename, const std::set< std::string > & columns, const std::vector< rowIndex_t > & rows, std::string delimiter = "\t", bool printRowNumbers = true)`

This function is equivalent in behavior to [print\(\)](#) except that instead of returning a string containing the printed contents of the data frame, the string is printed to a file. If the file exists, it will be overwritten.

Parameters

<i>filename</i>	Name of the file to print to. If it is an absolute path, the file will be put there. If it is a local path, the file will be placed relative to the data directory of the project.
-----------------	--

Returns

True for success, false if there was some problem writing to the file (insufficient permissions, etc.)

12.3.2.21 `bool CX_DataFrame::readFromFile (std::string filename, std::string cellDelimiter = "\t", std::string vectorEncloser = "\"\"\")`

Reads data from the given file into the data frame. This function assumes that there will be a row of column names as the first row of the file. It does not treat consecutive delimiters as a single delimiter.

Parameters

<i>filename</i>	The name of the file to read data from. If it is a relative path, the file will be read relative to the data directory.
<i>cellDelimiter</i>	A string containing the delimiter between cells of the data frame.
<i>vectorEncloser</i>	A string containing the characters that surround cell that contain a vector of data. By default, vectors are enclosed in double quotes. This indicates to most software that it should treat the contents of the quotes "as-is", i.e. if it finds a delimiter within the quotes, it should not split there, but wait until out of the quotes.

Returns

False if an error occurred, true otherwise.

Note

The contents of the data frame will be deleted before attempting to read in the file.
If the data is read in from a file written with a row numbers column, that column will be read into the data frame. You can remove it using `deleteColumn("rowNumber")`.

12.3.2.22 bool CX_DataFrame::reorderRows (const vector< CX_DataFrame::rowIndex_t > & newOrder)

Re-orders the rows in the data frame.

Parameters

<i>newOrder</i>	Vector of row indices. newOrder.size() must equal this->getRowCount(). newOrder must not contain any out-of-range indices (i.e. they must be < getRowCount()). Both of these error conditions are checked for in the function call and errors are logged.
-----------------	--

Returns

true if all of the conditions of newOrder are met, false otherwise.

12.3.2.23 void CX_DataFrame::setRowCount (rowIndex_t rowCount)

Sets the number of rows in the data frame.

Parameters

<i>rowCount</i>	The new number of rows in the data frame.
-----------------	---

Note

If the row count is less than the number of rows already in the data frame, it will delete those rows with a warning.

12.3.2.24 void CX_DataFrame::shuffleRows (void)

Randomly re-orders the rows of the data frame using [CX::Instances::RNG](#) as the random number generator for the shuffling.

12.3.2.25 void CX_DataFrame::shuffleRows (CX_RandomNumberGenerator & rng)

Randomly re-orders the rows of the data frame.

Parameters

<i>rng</i>	Reference to a CX_RandomNumberGenerator to be used for the shuffling.
------------	---

The documentation for this class was generated from the following files:

- CX_DataFrame.h
- CX_DataFrame.cpp

12.4 CX::CX_DataFrameCell Class Reference

```
#include <CX_DataFrameCell.h>
```

Public Member Functions

- [CX_DataFrameCell](#) (const char *c)
 - [CX_DataFrameCell](#) (const T &value)
- Construct the cell, assigning the value to it.*

- `template<typename T >`
`CX_DataFrameCell` (const std::vector< T > &values)
Construct the cell, assigning the values to it.
- `CX_DataFrameCell` & `operator=` (const char *c)
- `template<typename T >`
`CX_DataFrameCell` & `operator=` (const T &value)
Assigns a value to the cell.
- `template<typename T >`
`CX_DataFrameCell` & `operator=` (const std::vector< T > &values)
Assigns a vector of values to the cell.
- `template<typename T >`
`operator T` (void) const
Attempts to convert the contents of the cell to T using `to()`.
- `template<typename T >`
`operator std::vector< T >` (void) const
Attempts to convert the contents of the cell to vector<T> using `toVector<T>()`.
- `template<typename T >`
void `store` (const T &value)
- `template<typename T >`
T `to` (void) const
- `std::string` `toString` (void) const
Returns a copy of the stored data in the internal string representation. Type checking is not done because this is a lossless operation.
- bool `toBool` (void) const
Returns a copy of the stored data converted to bool. Equivalent to `to<bool>()`.
- int `toInt` (void) const
Returns a copy of the stored data converted to int. Equivalent to `to<int>()`.
- double `toDouble` (void) const
Returns a copy of the stored data converted to double. Equivalent to `to<double>()`.
- `template<typename T >`
std::vector< T > `toVector` (void) const
- `template<typename T >`
void `storeVector` (std::vector< T > values)
- void `copyCellTo` (CX_DataFrameCell *targetCell)
- std::string `getStoredType` (void)
- `template<>`
std::string `to` (void) const

12.4.1 Detailed Description

This class manages the contents of a single cell in a `CX_DataFrame`. It handles all of the type conversion nonsense that goes on when data is inserted into or extracted from a data frame. It tracks the type of the data that is inserted or extracted and logs warnings if the inserted type does not match the extracted type, with a few exceptions (see notes).

Note

There are a few exceptions to the type tracking. If the inserted type is `const char*`, it is treated as a string. Additionally, you can extract anything as string without a warning. This is because the data is stored as a string internally so extracting the data as a string is a lossless operation.

12.4.2 Constructor & Destructor Documentation

12.4.2.1 CX_DataFrameCell::CX_DataFrameCell (const char * c)

Constructs the cell with a string literal, treating it as a std::string.

12.4.3 Member Function Documentation

12.4.3.1 void CX_DataFrameCell::copyCellTo (CX_DataFrameCell * targetCell)

Copies the contents of this cell to targetCell, including type information.

Parameters

<i>targetCell</i>	A pointer to the cell to copy data to.
-------------------	--

12.4.3.2 std::string CX_DataFrameCell::getStoredType (void)

Gets a string representing the type of data stored within the cell. This string is implementation-defined (which is the C++ standards committee way of saying "It can be anything at all"). It is only guaranteed to be the same for the same type, but not necessarily be different for different types.

Returns

A string containing the name of the stored type as given by typeid(typename).name().

12.4.3.3 CX_DataFrameCell & CX_DataFrameCell::operator= (const char * c)

Assigns a string literal to the cell, treating it as a std::string.

12.4.3.4 template<typename T> void CX::CX_DataFrameCell::store (const T & value)

Stores the given value with the given type. This function is a good way to explicitly state the type of the data you are storing into the cell if, for example, it is a literal.

Template Parameters

<i>< T></i>	The type to store the value as. If T is not specified, this function is essentially equivalent to using operator=.
-------------------	--

Parameters

<i>value</i>	The value to store.
--------------	---------------------

12.4.3.5 template<typename T> void CX::CX_DataFrameCell::storeVector (std::vector< T> values)

Stores a vector of data in the cell. The data is stored as a string with each element delimited by a semicolon. If the data to be stored are strings containing semicolons, the data will not be extracted properly.

Parameters

<i>values</i>	A vector of values to store.
---------------	------------------------------

12.4.3.6 `template<typename T> T CX::CX_DataFrameCell::to (void) const`

Attempts to convert the contents of the cell to type T. There are a variety of reasons why this conversion can fail and they all center on the user inserting data of one type and then attempting to extract data of a different type. Regardless of whether the conversion is possible, if you try to extract a type that is different from the type that is stored in the cell, a warning will be logged.

Template Parameters

<code><T></code>	The type to convert to.
------------------------	-------------------------

Returns

The data in the cell converted to T.

12.4.3.7 `std::string CX::CX_DataFrameCell::to (void) const`

Equivalent to a call to [toString\(\)](#). This is specialized because it skips the type checks of `to<T>`.

Returns

A copy of the stored data encoded as a string.

12.4.3.8 `template<typename T> std::vector< T > CX::CX_DataFrameCell::toVector (void) const`

Returns a copy of the contents of the cell converted to a vector of the given type. If the type of data stored in the cell was not a vector of the given type or the type does match but it was a scalar that is stored, the logs a warning but attempts the conversion anyway.

Template Parameters

<code><T></code>	The type of the elements of the returned vector.
------------------------	--

Returns

A vector containing the converted data.

The documentation for this class was generated from the following files:

- CX_DataFrameCell.h
- CX_DataFrameCell.cpp

12.5 CX::CX_DataFrameColumn Class Reference

Public Member Functions

- [CX_DataFrameCell](#) **operator[]** (CX_DataFrame::RowIndex_t row)
- CX_DataFrame::RowIndex_t **size** (void)

Friends

- class **CX_DataFrame**

12.5.1 Detailed Description

The documentation for this class was generated from the following files:

- CX_DataFrame.h
- CX_DataFrame.cpp

12.6 CX::CX_DataFrameRow Class Reference

Public Member Functions

- [CX_DataFrameCell](#) **operator[]** (std::string column)
- vector< std::string > **names** (void)
- void **clear** (void)

Friends

- class **CX_DataFrame**

12.6.1 Detailed Description

The documentation for this class was generated from the following files:

- CX_DataFrame.h
- CX_DataFrame.cpp

12.7 CX::Util::CX_DegreeToPixelConverter Class Reference

```
#include <CX_UnitConversion.h>
```

Inherits CX::Util::CX_BaseUnitConverter.

Public Member Functions

- [CX_DegreeToPixelConverter](#) (float pixelsPerUnit, float viewingDistance, bool roundResult=false)
- float **operator()** (float degrees)

12.7.1 Detailed Description

This simple utility class is used for converting degrees of visual angle to pixels on a monitor. This class uses CX::Util::degreesToPixels() internally. See also [CX::Util::CX_CoordinateConverter](#) for a way to also convert from one coordinate system to another.

Example use:

```
CX_DegreeToPixelConverter d2p(34, 60); //34 pixels per unit length (e.g. cm) on
    the target monitor, user is 60 length units from monitor.
ofLine( 200, 100, 200 + d2p(1), 100 + d2p(2) ); //Draw a line from (200, 100) (in pixel coordinates) to 1
    degree
//to the right and 2 degrees below that point.
```

12.7.2 Constructor & Destructor Documentation

12.7.2.1 `CX::Util::CX_DegreeToPixelConverter::CX_DegreeToPixelConverter (float pixelsPerUnit, float viewingDistance, bool roundResult = false)`

Constructs an instance of a [CX_DegreeToPixelConverter](#) using the given settings.

Parameters

<i>pixelsPerUnit</i>	The number of pixels within one length unit (e.g. inches, centimeters). This can be measured by drawing a ~100-1000 pixel square on the screen and measuring the length of a side and dividing the number of pixels by the total length measured.
<i>viewingDistance</i>	The distance from the monitor that the participant will be viewing the screen from.
<i>roundResult</i>	If true, the result of conversions will be rounded to the nearest integer (i.e. pixel). For drawing certain kinds of stimuli (especially text) it can be helpful to draw on pixel boundaries.

12.7.3 Member Function Documentation

12.7.3.1 float CX::Util::CX_DegreeToPixelConverter::operator() (float *degrees*)

Converts the degrees to pixels based on the settings given during construction.

Parameters

<i>degrees</i>	The number of degrees of visual angle to convert to pixels.
----------------	---

Returns

The number of pixels corresponding to the number of degrees of visual angle.

The documentation for this class was generated from the following files:

- CX_UnitConversion.h
- CX_UnitConversion.cpp

12.8 CX::CX_Display Class Reference

```
#include <CX_Display.h>
```

Public Member Functions

- void [setup](#) (void)
- void [setFullScreen](#) (bool fullScreen)
- void [drawFboToBackBuffer](#) (ofFbo &fbo)
- void [drawFboToBackBuffer](#) (ofFbo &fbo, ofRectangle placement)
- void [beginDrawingToBackBuffer](#) (void)
- void [endDrawingToBackBuffer](#) (void)
- void [BLOCKING_swapFrontAndBackBuffers](#) (void)
- void [swapFrontAndBackBuffers](#) (void)
- void [BLOCKING_setAutoSwapping](#) (bool autoSwap)
- bool [isAutomaticallySwapping](#) (void)
- bool [hasSwappedSinceLastCheck](#) (void)
- CX_Micros [getLastSwapTime](#) (void)
- CX_Micros [getFramePeriod](#) (void)
- void [setWindowResolution](#) (int width, int height)
- void [setWindowTitle](#) (std::string title)
- ofRectangle [getResolution](#) (void)
- ofPoint [getCenterOfDisplay](#) (void)

- `uint64_t getFrameNumber` (void)
- `void BLOCKING_estimateFramePeriod` (CX_Micros estimationInterval)
- `CX_Micros estimateNextSwapTime` (void)
- `void BLOCKING_waitForOpenGL` (void)

12.8.1 Detailed Description

This class represents an abstract visual display surface, which is my way of saying that it doesn't necessarily represent a monitor. The display surface can either be a window or, if full screen, the whole monitor. It is also a bit abstract in that it does not draw anything, but only creates an context in which things can be drawn.

12.8.2 Member Function Documentation

12.8.2.1 `void CX_Display::beginDrawingToBackBuffer (void)`

Prepares a rendering context for using drawing functions. Must be paired with a call to `endDrawingToBackBuffer()`.

12.8.2.2 `void CX_Display::BLOCKING_estimateFramePeriod (CX_Micros estimationInterval)`

This function estimates the typical period of the display refresh. This function blocks for estimationInterval while the swapping thread swaps in the background. This function is called with an argument of 300 ms during construction of this class, so there will always be some information about the frame period. If more precision of the estimate is desired, this function can be called again with a longer wait duration.

Parameters

<i>estimationInterval</i>	The length of time to spend estimating the frame period.
---------------------------	--

12.8.2.3 `void CX_Display::BLOCKING_setAutoSwapping (bool autoSwap)`

Set whether the front and buffers of the display will swap automatically every frame or not. You can check to see if a swap has occurred by calling `hasSwappedSinceLastCheck()`. You can check to see if the display is automatically swapping by calling `isAutomaticallySwapping()`.

Parameters

<i>autoSwap</i>	If true, the front and back buffer will swap automatically every frame.
-----------------	---

12.8.2.4 `void CX_Display::BLOCKING_swapFrontAndBackBuffers (void)`

This function queues up a swap of the front and back buffers then blocks until the swap occurs. It does nothing if `isAutomaticallySwapping() == true`.

12.8.2.5 `void CX_Display::BLOCKING_waitForOpenGL (void)`

Wait until all OpenGL instructions that were given before this was called to complete. Any commands put into the pipeline after this is called (from other threads) are not waited for.

12.8.2.6 `void CX_Display::drawFboToBackBuffer (ofFbo & fbo)`

Draw the given ofFbo to the back buffer. It will be drawn starting from 0, 0 and will be drawn at the full dimensions of the ofFbo (whatever size was chosen at allocation of the fbo).

12.8.2.7 void CX_Display::drawFboToBackBuffer (ofFbo & *fbo*, ofRectangle *rect*)

[Draw](#) the given ofFbo to the back buffer at the coordinates given by rect.

Parameters

<i>fbo</i>	The fbo to draw.
<i>rect</i>	The rectangle in which to place to fbo. The x and y components specify location. The width and height components specify the output width and height of the fbo. If these are not equal to the width and height of the fbo, the fbo will be scaled up or down to fit the width and height.

12.8.2.8 void CX_Display::endDrawingToBackBuffer (void)

Finish rendering to the back buffer. Must be paired with a call to [beginDrawingToBackBuffer\(\)](#).

12.8.2.9 CX_Micros CX_Display::estimateNextSwapTime (void)

Get an estimate of the next time the front and back buffers will be swapped. This function depends on the precision of the frame period as estimated using [BLOCKING_estimateFramePeriod\(\)](#).

Returns

A time value that can be compared to CX::Instances::Clock.getTime().

12.8.2.10 ofPoint CX_Display::getCenterOfDisplay (void)

Returns an ofPoint representing the center of the display. Works in either windowed or full screen mode.

12.8.2.11 uint64_t CX_Display::getFrameNumber (void)

This function returns the number of the last frame presented, as determined by number of front and back buffer swaps. It tracks buffer swaps that result from 1) the front and back buffer swapping automatically (as a result of [BLOCKING_setAutoSwapping\(\)](#) with true as the argument) and 2) manual swaps resulting from a call to [BLOCKING_swapFrontAndBackBuffers\(\)](#) or [swapFrontAndBackBuffers\(\)](#).

Returns

The number of the last frame. This value can only be compared with other values returned by this function.

12.8.2.12 CX_Micros CX_Display::getFramePeriod (void)

Gets the estimate of the frame period calculated with [BLOCKING_estimateFramePeriod\(\)](#).

12.8.2.13 CX_Micros CX_Display::getLastSwapTime (void)

Get the last time at which the front and back buffers were swapped.

Returns

A time value that can be compared with CX::Instances::Clock.getTime().

12.8.2.14 ofRectangle CX_Display::getResolution (void)

Returns the resolution of the current window, not the resolution of the monitor (unless you are in full screen mode). You can use either x and y or width and height.

12.8.2.15 bool CX_Display::hasSwappedSinceLastCheck (void)

Check to see if the display has swapped the front and back buffers since the last call to this function. This is generally used in conjunction with automatic swapping of the buffers ([BLOCKING_setAutoSwapping\(\)](#)) or with an individual

threaded swap of the buffers ([swapFrontAndBackBuffers\(\)](#)). This technically works with [BLOCKING_swapFrontAndBackBuffers\(\)](#), but given that that function only returns once the buffers have swapped, checking that the buffers have swapped is redundant.

Returns

True if a swap has been made since the last call to this function, false otherwise.

12.8.2.16 bool CX_Display::isAutomaticallySwapping (void)

Determine whether the display is configured to automatically swap the front and back buffers every frame. See [BLOCKING_setAutoSwapping](#) for more information.

12.8.2.17 void CX_Display::setFullScreen (bool *fullScreen*)

Set whether the display is full screen or not. If the display is set to full screen, the resolution may not be the same as the resolution of display in windowed mode, and vice versa.

12.8.2.18 void CX_Display::setup (void)

Set up the display. Must be called for the display to function correctly.

12.8.2.19 void CX_Display::setWindowResolution (int *width*, int *height*)

Sets the resolution of the window. Has no effect if called while in full screen mode.

Parameters

<i>width</i>	The desired width of the window, in pixels.
<i>height</i>	The desired height of the window, in pixels.

12.8.2.20 void CX_Display::setWindowTitle (std::string *title*)

Sets the title of the experiment window.

Parameters

<i>title</i>	The new window title.
--------------	-----------------------

12.8.2.21 void CX_Display::swapFrontAndBackBuffers (void)

This function cues a swap of the front and back buffers. It avoids blocking (like [BLOCKING_swapFrontAndBackBuffers\(\)](#)) by spawning a thread in which the swap is waited for.

The documentation for this class was generated from the following files:

- CX_Display.h
- CX_Display.cpp

12.9 CX::CX_FinalSlideFunctionInfo_t Struct Reference

```
#include <CX_SlidePresenter.h>
```

Public Attributes

- [CX_SlidePresenter](#) * *instance*

A pointer to the [CX_SlidePresenter](#) that called the user function.

- unsigned int [currentSlideIndex](#)

The index of the slide that is currently being presented.

12.9.1 Detailed Description

The final slide function takes a reference to a struct of this type.

The documentation for this struct was generated from the following file:

- [CX_SlidePresenter.h](#)

12.10 CX::CX_InputManager Class Reference

```
#include <CX_InputManager.h>
```

Public Member Functions

- bool [setup](#) (bool useKeyboard, bool useMouse, int joystickIndex=-1)
- bool [pollEvents](#) (void)

Public Attributes

- [CX_Keyboard](#) Keyboard
An instance of [CX::CX_Keyboard](#). Configured using [CX::CX_InputManager::setup\(\)](#).
- [CX_Mouse](#) Mouse
An instance of [CX::CX_Mouse](#). Configured using [CX::CX_InputManager::setup\(\)](#).
- [CX_Joystick](#) Joystick
An instance of [CX::CX_Joystick](#). Configured using [CX::CX_InputManager::setup\(\)](#).

12.10.1 Detailed Description

This class is responsible for managing three basic input devices: a keyboard, mouse, and joystick. You access each of these devices with the corresponding member class: Keyboard, Mouse, and Joystick. See [CX::CX_Keyboard](#), [CX::CX_Mouse](#), and [CX::CX_Joystick](#) for more information about each specific device.

By default, all three input devices are disabled. Call [setup\(\)](#) to enable specific devices.

12.10.2 Member Function Documentation

12.10.2.1 bool CX_InputManager::pollEvents (void)

It is not typically necessary for the user to call this function directly, although there is no harm in doing so. This function polls for new events on all of the configured input devices (see [setup\(\)](#)). After a call to this function, new events for the input devices can be found by checking the [availableEvents\(\)](#) function for each device.

Returns

True if there are any new events, false otherwise.

12.10.2.2 bool CX_InputManager::setup (bool *useKeyboard*, bool *useMouse*, int *joystickIndex* = -1)

Setup the input manager to use the requested devices. You may call this function multiple times if you want to change the configuration over the course of the experiment. Every time this function is called, all input device events are cleared.

Parameters

<i>useKeyboard</i>	Enable or disable the keyboard.
<i>useMouse</i>	Enable or disable the mouse.
<i>joystickIndex</i>	Optional. If ≥ 0 , an attempt will be made to set up the joystick at that index. If < 0 , no attempt will be made to set up the joystick.

Returns

False if the requested joystick could not be set up correctly, true otherwise.

The documentation for this class was generated from the following files:

- CX_InputManager.h
- CX_InputManager.cpp

12.11 CX::CX_Joystick Class Reference

```
#include <CX_Joystick.h>
```

Public Member Functions

- bool [setup](#) (int joystickIndex)
- std::string [getJoystickName](#) (void)
- bool [pollEvents](#) (void)
- int [availableEvents](#) (void)
- CX_JoystickEvent_t [getNextEvent](#) (void)
- void [clearEvents](#) (void)
- std::vector< float > [getAxisPositions](#) (void)
- std::vector< unsigned char > [getButtonStates](#) (void)

12.11.1 Detailed Description

This class manages a joystick that is attached to the system (if any). If more than one joystick is needed for the experiment, you can create more instances of [CX_Joystick](#) other than the one in [CX::Instances::Input](#).

12.11.2 Member Function Documentation**12.11.2.1 int CX_Joystick::availableEvents (void)**

Get the number of new events available for this input device.

12.11.2.2 void CX_Joystick::clearEvents (void)

Clear (delete) all events from this input device.

12.11.2.3 vector< float > CX_Joystick::getAxisPositions (void)

This function is to be used for direct access to the axis positions of the joystick. It does not generate events (i.e. [CX_JoystickEvent_t](#)), nor does it do any timestamping. If timestamps and uncertainties are desired, you MUST use [pollEvents\(\)](#) and the associated event functions (e.g. [getNextEvent\(\)](#)).

12.11.2.4 `vector< unsigned char > CX_Joystick::getButtonStates (void)`

This function is to be used for direct access to the button states of the joystick. It does not generate events (i.e. [CX_JoystickEvent_t](#)), nor does it do any timestamping. If timestamps and uncertainties are desired, you MUST use [pollEvents\(\)](#) and the associated event functions (e.g. [getNextEvent\(\)](#)).

12.11.2.5 `std::string CX_Joystick::getJoystickName (void)`

Get the name of the joystick, presumably as set by the joystick driver. The name may not be very meaningful.

12.11.2.6 `CX_JoystickEvent_t CX_Joystick::getNextEvent (void)`

Get the next event available for this input device. This is a destructive operation: the returned event is deleted from the input device.

12.11.2.7 `bool CX_Joystick::pollEvents (void)`

Check to see if there are any new joystick events. If there are new events, they can be accessed with [availableEvents\(\)](#) and [getNextEvent\(\)](#).

Returns

True if there are new events.

12.11.2.8 `bool CX_Joystick::setup (int joystickIndex)`

Set up the joystick by attempting to initialize the joystick at the given index. If the joystick is present on the system, it will be initialized and its name can be accessed by calling [getJoystickName\(\)](#).

If the set up is successful (i.e. if the selected joystick is present on the system), this function will return true. If the joystick is not present, it will return false.

The documentation for this class was generated from the following files:

- CX_Joystick.h
- CX_Joystick.cpp

12.12 CX::CX_JoystickEvent_t Struct Reference

```
#include <CX_Joystick.h>
```

Public Types

- enum [JoystickEventType](#) { [BUTTON_PRESS](#), [BUTTON_RELEASE](#), [AXIS_POSITION_CHANGE](#) }

Public Attributes

- int [buttonIndex](#)
If eventType is [BUTTON_PRESS](#) or [BUTTON_RELEASE](#), this contains the index of the button that was changed.
- unsigned char [buttonState](#)
If eventType is [BUTTON_PRESS](#) or [BUTTON_RELEASE](#), this contains the current state of the button.
- int [axisIndex](#)
If eventType is [AXIS_POSITION_CHANGE](#), this contains the index of the axis which changed.

- float [axisPosition](#)
If eventType is `AXIS_POSITION_CHANGE`, this contains the amount by which the axis changed.
- CX_Micros [eventTime](#)
The time at which the event was registered. Can be compared to the result of `CX::Clock::getTime()`.
- CX_Micros [uncertainty](#)
The uncertainty in eventTime. The event occurred some time between eventTime and eventTime minus uncertainty.
- enum
[CX::CX_JoystickEvent_t::JoystickEventType eventType](#)
The type of the event.

12.12.1 Detailed Description

This struct contains information about joystick events. Joystick events are either a button press or release or a change in the axes of the joystick.

See Also

[CX::CX_Joystick::getNextEvent\(\)](#) provides access to joystick events.

12.12.2 Member Enumeration Documentation

12.12.2.1 enum CX::CX_JoystickEvent_t::JoystickEventType

Enumerator

- BUTTON_PRESS** A button on the joystick has been pressed. See [buttonIndex](#) and [buttonState](#) for the event data.
- BUTTON_RELEASE** A button on the joystick has been released. See [buttonIndex](#) and [buttonState](#) for the event data.
- AXIS_POSITION_CHANGE** The joystick has been moved in one of its axes. See [axisIndex](#) and [axisPosition](#) for the event data.

The documentation for this struct was generated from the following file:

- CX_Joystick.h

12.13 CX::CX_Keyboard Class Reference

```
#include <CX_Keyboard.h>
```

Public Member Functions

- int [availableEvents](#) (void)
- [CX_KeyEvent_t getNextEvent](#) (void)
- void [clearEvents](#) (void)

Friends

- class [CX_InputManager](#)

12.13.1 Detailed Description

This class is responsible for managing the mouse.

12.13.2 Member Function Documentation

12.13.2.1 int CX_Keyboard::availableEvents (void)

Get the number of new events available for this input device.

12.13.2.2 void CX_Keyboard::clearEvents (void)

Clear (delete) all events from this input device.

12.13.2.3 CX_KeyEvent_t CX_Keyboard::getNextEvent (void)

Get the next event available for this input device. This is a destructive operation: the returned event is deleted from the input device.

The documentation for this class was generated from the following files:

- CX_Keyboard.h
- CX_Keyboard.cpp

12.14 CX::CX_KeyEvent_t Struct Reference

```
#include <CX_Keyboard.h>
```

Public Types

- enum [KeyboardEventType](#) { [PRESSED](#), [RELEASED](#), [REPEAT](#) }

Public Attributes

- int [key](#)
- CX_Micros [eventTime](#)
The time at which the event was registered. Can be compared to the result of CX::Clock::getTime().
- CX_Micros [uncertainty](#)
The uncertainty in eventTime. The event occurred some time between eventTime and eventTime minus uncertainty.
- enum [CX::CX_KeyEvent_t::KeyboardEventType](#) [eventType](#)
The type of the event.

12.14.1 Detailed Description

This struct contains the results of a keyboard event, whether it be a key press or release, or key repeat.

12.14.2 Member Enumeration Documentation

12.14.2.1 enum CX::CX_KeyEvent_t::KeyboardEventType

Enumerator

PRESSED A key has been pressed.

RELEASED A key has been released.

REPEAT A key has been held for some time and automatic key repeat has kicked in, causing multiple keypresses to be rapidly sent.

12.14.3 Member Data Documentation

12.14.3.1 int CX::CX_KeyEvent_t::key

The key involved in this event. The value of this can be compared with chars for many keys (e.g. (myKeyEvent.key == 'e') to test if the key was the E key. For special keys, this can be compared with values defined in ofConstants.h (e.g. OF_KEY_ESC).

The documentation for this struct was generated from the following file:

- CX_Keyboard.h

12.15 CX::Util::CX_LengthToPixelConverter Class Reference

```
#include <CX_UnitConversion.h>
```

Inherits CX::Util::CX_BaseUnitConverter.

Public Member Functions

- [CX_LengthToPixelConverter](#) (float pixelsPerUnit, bool roundResult=false)
- float [operator\(\)](#) (float length)

12.15.1 Detailed Description

This simple utility class is used for converting lengths (perhaps of objects drawn on the monitor) to pixels on a monitor. See also [CX::Util::CX_CoordinateConverter](#) for a way to also convert from one coordinate system to another.

Example use:

```
CX_LengthToPixelConverter l2p(75); //75 pixels per unit length (e.g. inch) on the
    target monitor.
ofLine( 200, 100, 200 + l2p(1), 100 + l2p(2) ); //Draw a line from (200, 100) (in pixel coordinates) to 1
    unit
//horizontally and 2 units vertically from that point.
```

12.15.2 Constructor & Destructor Documentation

12.15.2.1 CX::Util::CX_LengthToPixelConverter::CX_LengthToPixelConverter (float pixelsPerUnit, bool roundResult = false)

Constructs a [CX_LengthToPixelConverter](#) with the given configuration.

Parameters

<i>pixelsPerUnit</i>	The number of pixels per one length unit. This can be measured by drawing a ~100-1000 pixel square on the screen and measuring the length of a side and dividing the number of pixels by the total length measured.
<i>roundResult</i>	If true, the result of conversions will be rounded to the nearest integer (i.e. pixel). For drawing certain kinds of stimuli (especially text) it can be helpful to draw on pixel boundaries.

12.15.3 Member Function Documentation

12.15.3.1 float CX::Util::CX_LengthToPixelConverter::operator() (float *length*)

Converts the length to pixels based on the settings given during construction.

Parameters

<i>length</i>	The length to convert to pixels.
---------------	----------------------------------

Returns

The number of pixels corresponding to the length.

The documentation for this class was generated from the following files:

- CX_UnitConversion.h
- CX_UnitConversion.cpp

12.16 CX::CX_Logger Class Reference

```
#include <CX_Logger.h>
```

Public Member Functions

- `std::stringstream & log (CX_LogLevel level, std::string module="")`
- `std::stringstream & verbose (std::string module="")`
- `std::stringstream & notice (std::string module="")`
- `std::stringstream & warning (std::string module="")`
- `std::stringstream & error (std::string module="")`
- `std::stringstream & fatalError (std::string module="")`
- `void level (CX_LogLevel level, std::string module="")`
- `void levelForConsole (CX_LogLevel level)`
- `void levelForFile (CX_LogLevel level, std::string filename="CX_DEFERRED_LOGGER_DEFAULT")`
- `void levelForAllModules (CX_LogLevel level)`
- `void flush (void)`
- `void timestamps (bool logTimestamps, std::string format="%H:%M:%S.%i")`
- `void setMessageFlushCallback (std::function< void(CX_MessageFlushData &)> f)`
- `void captureOFLogMessages (void)`

12.16.1 Detailed Description

This class is used for logging messages throughout the [CX](#) backend code. It can also be used in user code to log messages. Rather than instantiating your own copy of [CX_Logger](#), it is probably better to use the preinstantiated [CX::Instances::Log](#).

There is an example showing a number of the features of [CX_Logger](#) named example-logging.

12.16.2 Member Function Documentation

12.16.2.1 void CX_Logger::captureOFLogMessages (void)

Set this instance of [CX_Logger](#) to be the target of any messages created by oF logging functions.

12.16.2.2 std::stringstream & CX_Logger::error (std::string module = " ")

This function is equivalent to a call to `log(CX_LogLevel::LOG_ERROR, module)`.

12.16.2.3 std::stringstream & CX_Logger::fatalError (std::string module = " ")

This function is equivalent to a call to `log(CX_LogLevel::LOG_FATAL_ERROR, module)`.

12.16.2.4 void CX_Logger::flush (void)

Log all of the messages stored since the last call to [flush\(\)](#) to the selected logging targets. This is a BLOCKING operation.

12.16.2.5 void CX_Logger::level (CX_LogLevel level, std::string module = " ")

Sets the log level for the given module. Messages from that module that are at a lower level than [level](#) will be ignored.

Parameters

<i>level</i>	See the CX::CX_LogLevel enum for valid values.
<i>module</i>	A string representing one of the modules from which log messages are generated.

12.16.2.6 void CX_Logger::levelForAllModules (CX_LogLevel level)

Set the log level for all modules. This works both retroactively and proactively: All currently known modules are given the log level and the default log level for new modules as set to the level.

12.16.2.7 void CX_Logger::levelForConsole (CX_LogLevel level)

Set the log level for messages to be printed to the console.

12.16.2.8 void CX_Logger::levelForFile (CX_LogLevel level, std::string filename = "CX_DEFERRED_LOGGER_DEFAULT")

Sets the log level for the file with given file name. If the file does not exist, it will be created. If the file does exist, it will be overwritten with a warning logged to cerr.

Parameters

<i>level</i>	See the CX_LogLevel enum for valid values.
<i>filename</i>	Optional. If no file name is given, a file with name generated from a date/time from the start time of the experiment will be used.

12.16.2.9 std::stringstream & CX_Logger::log (CX_LogLevel *level*, std::string *module* = " ")

This is the basic logging function for this class. Example use:

```
Log.log(CX_LogLevel::LOG_WARNING, "myModule") << "My message number " << 20;
```

Possible output: "[warning] <myModule> My message number 20"

Parameters

<i>level</i>	Log level for this message. This has implications for message filtering. See level() . This should not be LOG_ALL or LOG_NONE, because that would be weird, wouldn't it?
<i>module</i>	Name of the module that this log message is related to. This has implications for message filtering. See level() .

Returns

A reference to a std::stringstream that the log message data should be streamed into.

12.16.2.10 std::stringstream & CX_Logger::notice (std::string *module* = " ")

This function is equivalent to a call to log(CX_LogLevel::LOG_NOTICE, module).

12.16.2.11 void CX_Logger::setMessageFlushCallback (std::function< void(CX_MessageFlushData &)> *f*)

Sets the user function that will be called on each message flush event. For every message that has been logged, the user function will be called. No filtering is performed: All messages regardless of the module log level will be sent to the user function.

Parameters

<i>f</i>	A pointer to a user function that takes a reference to a CX_MessageFlushData struct and returns nothing.
----------	--

12.16.2.12 void CX_Logger::timestamps (bool *logTimestamps*, std::string *format* = "%H:%M:%S.%i")

Set whether or not to log timestamps and the format for the timestamps.

Parameters

<i>logTimestamps</i>	Does what it says.
<i>format</i>	Timestamp format string. See http://pocoproject.org/docs/Poco.DateTime-Formatter.html#4684 for documentation of the format. Defaults to H:M:S.i (24-hour clock with milliseconds at the end).

12.16.2.13 std::stringstream & CX_Logger::verbose (std::string *module* = " ")

This function is equivalent to a call to log(CX_LogLevel::LOG_VERBOSE, module).

12.16.2.14 std::stringstream & CX_Logger::warning (std::string *module* = " ")

This function is equivalent to a call to log(CX_LogLevel::LOG_WARNING, module).

The documentation for this class was generated from the following files:

- CX_Logger.h
- CX_Logger.cpp

12.17 CX::CX_Mouse Class Reference

```
#include <CX_Mouse.h>
```

Public Member Functions

- int [availableEvents](#) (void)
- [CX_MouseEvent_t getNextEvent](#) (void)
- void [clearEvents](#) (void)
- void [showCursor](#) (bool show)
- void [setCursorPosition](#) (ofPoint pos)
- ofPoint [getCursorPosition](#) (void)

Friends

- class **CX_InputManager**

12.17.1 Detailed Description

This class is responsible for managing the mouse.

12.17.2 Member Function Documentation

12.17.2.1 int CX_Mouse::availableEvents (void)

Get the number of new events available for this input device.

12.17.2.2 void CX_Mouse::clearEvents (void)

Clear (delete) all events from this input device.

12.17.2.3 ofPoint CX_Mouse::getCursorPosition (void)

Get the cursor position within the program window. If the mouse has left the window, this will return the last known position of the cursor within the window.

Returns

An ofPoint with the last cursor position.

12.17.2.4 CX_MouseEvent_t CX_Mouse::getNextEvent (void)

Get the next event available for this input device. This is a destructive operation: the returned event is deleted from the input device.

12.17.2.5 void CX_Mouse::setCursorPosition (ofPoint pos)

Sets the position of the cursor, relative to the program the window. The window must be focused.

Parameters

<i>pos</i>	The location within the window to set the cursor.
------------	---

12.17.2.6 void CX_Mouse::showCursor (bool show)

Show or hide the mouse cursor within the program window. If in windowed mode, the cursor will be visible outside of the window.

Parameters

<i>show</i>	If true, the cursor will be shown, if false it will not be shown.
-------------	---

The documentation for this class was generated from the following files:

- CX_Mouse.h
- CX_Mouse.cpp

12.18 CX::CX_MouseEvent_t Struct Reference

```
#include <CX_Mouse.h>
```

Public Types

- enum [MouseEventType](#) {
[MOVED](#), [PRESSED](#), [RELEASED](#), [DRAGGED](#),
[SCROLLED](#) }

Public Attributes

- int [button](#)
The relevant mouse button if the eventType is PRESSED, RELEASED, or DRAGGED. Can be compared with elements of enum CX_MouseButtons to find out about the primary buttons.
- int [x](#)
The x position of the cursor at the time of the event, or the change in the x-axis scroll if the eventType is SCROLLED.
- int [y](#)
The y position of the cursor at the time of the event, or the change in the y-axis scroll if the eventType is SCROLLED.
- CX_Micros [eventTime](#)
The time at which the event was registered. Can be compared to the result of CX::Clock::getTime().
- CX_Micros [uncertainty](#)
The uncertainty in eventTime. The event occurred some time between eventTime and eventTime minus uncertainty.
- enum
[CX::CX_MouseEvent_t::MouseEventType eventType](#)
The type of the event.

12.18.1 Detailed Description

This struct contains the results of a mouse event, which is any type of interaction with the mouse, be it simply movement, a button press or release, a drag event (mouse button held while mouse is moved), or movement of the scroll wheel.

12.18.2 Member Enumeration Documentation

12.18.2.1 enum CX::CX_MouseEvent_t::MouseEventType

Enumerator

MOVED The mouse has been moved without a button being held. [button](#) should be -1 (meaningless).

PRESSED A mouse button has been pressed. Check [button](#) for the button index and [x](#) and [y](#) for the location.

RELEASED A mouse button has been released. Check [button](#) for the button index and [x](#) and [y](#) for the location.

DRAGGED The mouse has been moved while at least one button was held. [button](#) may not be meaningful because the held button can be changed during a drag, or multiple buttons may be held at once during a drag.

SCROLLED The mouse wheel has been scrolled. Check [y](#) to get the change in the standard mouse wheel, or [x](#) if your mouse has a wheel that can move horizontally.

The documentation for this struct was generated from the following file:

- CX_Mouse.h

12.19 CX::CX_RandomNumberGenerator Class Reference

```
#include <CX_RandomNumberGenerator.h>
```

Public Member Functions

- [CX_RandomNumberGenerator](#) (void)
- void [setSeed](#) (unsigned long seed)
- unsigned long [getSeed](#) (void)
- CX_RandomInt_t [getMinimumRandomInt](#) (void)
- CX_RandomInt_t [getMaximumRandomInt](#) (void)
- CX_RandomInt_t [randomInt](#) (void)
- CX_RandomInt_t [randomInt](#) (CX_RandomInt_t rangeLower, CX_RandomInt_t rangeUpper)
- template<typename T >
T [randomExclusive](#) (const std::vector< T > &values, const T &exclude)
- template<typename T >
T [randomExclusive](#) (const std::vector< T > &values, const std::vector< T > &exclude)
- double [uniformDeviate](#) (double lowerBound_closed, double upperBound_open)
- std::vector< double > [uniformDeviates](#) (unsigned int count, double lowerBound_closed, double upperBound_open)
- template<typename T >
std::vector< T > [binomialDeviates](#) (unsigned int count, T trials, double probSuccess)
- std::vector< double > [normalDeviates](#) (unsigned int count, double mean, double standardDeviation)
- template<typename T >
void [shuffleVector](#) (std::vector< T > *v)
- template<typename T >
std::vector< T > [shuffleVector](#) (std::vector< T > v)
- template<typename T >
T [sample](#) (std::vector< T > values)
- template<typename T >
std::vector< T > [sample](#) (unsigned int count, const std::vector< T > &source, bool withReplacement)
- std::vector< int > [sample](#) (unsigned int count, int lowerBound, int upperBound, bool withReplacement)

12.19.1 Detailed Description

This class is used for generating random values from a pseudo-random number generator. It uses a version of the Mersenne Twister algorithm, in particular `std::mt19937_64` (see http://en.cppreference.com/w/cpp/numeric/random/mersenne_twister_engine for the parameters used with this algorithm).

The monolithic structure of `CX_RandomNumberGenerator` provides a certain important feature that a collection of loose function does not have, which is the ability to trivially track the random seed being used for the random number generator. The function `CX_RandomNumberGenerator::setSeed()` sets the seed for all random number generation tasks performed by this class. Likewise, `CX_RandomNumberGenerator::getSeed()` allows you to easily find the seed that is being used for random number generation. Due to this structure, you can easily save the seed that was used for each participant, which allows you to repeat the exact randomizations used for that participant (unless random number generation varies as a function of the responses given by a participant).

An instance of this class is preinstantiated for you. See `CX::Instances::RNG` for information about the instance.

12.19.2 Constructor & Destructor Documentation

12.19.2.1 CX_RandomNumberGenerator::CX_RandomNumberGenerator (void)

Constructs an instance of a `CX_RandomNumberGenerator`. Seeds the `CX_RandomNumberGenerator` using a `std::random_device`.

By the C++11 specification, `std::random_device` is supposed to be a non-deterministic (hardware) RNG. However, from http://en.cppreference.com/w/cpp/numeric/random/random_device: "Note that `std::random_device` may be implemented in terms of a pseudo-random number engine if a non-deterministic source (e.g. a hardware device) is not available to the implementation." According to a Stack Overflow comment, Microsoft's implementation of `std::random_device` is based on a ton of stuff, which should result in a fairly random result to be used as a seed for our Mersenne Twister. See the comment: <http://stackoverflow.com/questions/9549357/the-implementation-of-random-device-in-vs2010/9575747#9575747> Although this data should have high entropy, it is not a hardware RNG. The `random_device` is only used to seed the Mersenne Twister, so as long as the initial value is random enough, it should be fine.

12.19.3 Member Function Documentation

12.19.3.1 template<typename T> std::vector< T> CX::CX_RandomNumberGenerator::binomialDeviates (unsigned int count, T trials, double probSuccess)

Samples count deviates from a binomial distribution with the given number of trials and probability of success on each trial.

Parameters

<i>count</i>	The number of deviates to generate.
<i>trials</i>	The number of trials. Must be a non-negative integer.
<i>probSuccess</i>	The probability of a success on a given trial, where a success is the value 1.

Returns

A vector of the deviates.

12.19.3.2 CX_RandomInt_t CX_RandomNumberGenerator::getMaximumRandomInt (void)

Get the maximum possible value that can be returned by `randomInt()`.

Returns

The maximum value.

12.19.3.3 CX_RandomInt_t CX_RandomNumberGenerator::getMinimumRandomInt (void)

Get the minimum value that can be returned by [randomInt\(\)](#).

Returns

The minimum value.

12.19.3.4 unsigned long CX_RandomNumberGenerator::getSeed (void)

Get the seed used to seed the random number generator.

Returns

The seed. May have been set by the user with [setSeed\(\)](#) or during construction of the [CX_RandomNumber-Generator](#).

12.19.3.5 std::vector< double > CX_RandomNumberGenerator::normalDeviates (unsigned int *count*, double *mean*, double *standardDeviation*)

Samples count deviates from a normal distribution with the given mean and standard deviation.

Parameters

<i>count</i>	The number of deviates to generate.
<i>mean</i>	The mean of the distribution.
<i>standard-Deviation</i>	The standard deviation of the distribution.

Returns

A vector of the deviates.

12.19.3.6 template<typename T > T CX::CX_RandomNumberGenerator::randomExclusive (const std::vector< T > & *values*, const T & *exclude*)

Get a random value from a vector, without the possibility of getting the excluded value.

Parameters

<i>values</i>	The vectors of values to sample from.
<i>exclude</i>	The value to exclude from sampling.

Returns

The sampled value.

Note

If all of the values are excluded, an error will be logged and T() will be returned.

12.19.3.7 `template<typename T> T CX::CX_RandomNumberGenerator::randomExclusive (const std::vector< T > & values,
const std::vector< T > & exclude)`

Get a random value from a vector without the possibility of getting any of the excluded values.

Parameters

<i>values</i>	The vector of values to sample from.
<i>exclude</i>	The vector of values to exclude from sampling.

Returns

The sampled value.

Note

If all of the values are excluded, an error will be logged and `T()` will be returned.

12.19.3.8 `CX_RandomInt_t CX_RandomNumberGenerator::randomInt (void)`

Get a random integer in the range `getMinimumRandomInt()`, `getMaximumRandomInt()`, inclusive.

Returns

The int.

12.19.3.9 `CX_RandomInt_t CX_RandomNumberGenerator::randomInt (CX_RandomInt_t min, CX_RandomInt_t max)`

This function returns an integer from the range `[rangeLower, rangeUpper]`. The minimum and maximum values for the int returned from this function are given by `getMinimumRandomInt()` and `getMaximumRandomInt()`.

If `rangeLower > rangeUpper`, the lower and upper ranges are swapped. If `rangeLower == rangeUpper`, it returns `rangeLower`.

12.19.3.10 `template<typename T> T CX::CX_RandomNumberGenerator::sample (std::vector< T > values)`

Returns a single value sampled randomly from values.

Returns

The sampled value.

Note

If `values.size() == 0`, an error will be logged and `T()` will be returned.

12.19.3.11 `template<typename T> std::vector< T > CX::CX_RandomNumberGenerator::sample (unsigned int count, const std::vector< T > & source, bool withReplacement)`

Returns a vector of count values drawn randomly from source, with or without replacement. The returned values are in a random order.

Parameters

<i>count</i>	The number of samples to draw.
--------------	--------------------------------

<i>source</i>	A vector to be sampled from.
<i>withReplacement</i>	Sample with or without replacement.

Returns

A vector of the sampled values.

Note

If (count > source.size() && withReplacement == false), an empty vector is returned.

12.19.3.12 `std::vector< int > CX_RandomNumberGenerator::sample (unsigned int count, int lowerBound, int upperBound, bool withReplacement)`

Returns a vector of count integers drawn randomly from the range [lowerBound, upperBound] with or without replacement.

Parameters

<i>count</i>	The number of samples to draw.
<i>lowerBound</i>	The lower bound of the range to sample from. It is possible to sample this value.
<i>upperBound</i>	The upper bound of the range to sample from. It is possible to sample this value.
<i>withReplacement</i>	Sample with or without replacement.

Returns

A vector of the samples.

12.19.3.13 `void CX_RandomNumberGenerator::setSeed (unsigned long seed)`

Set the seed for the random number generator. You can retrieve the seed with [getSeed\(\)](#).

Parameters

<i>seed</i>	The new seed.
-------------	---------------

12.19.3.14 `template<typename T> void CX::CX_RandomNumberGenerator::shuffleVector (std::vector< T > * v)`

Randomizes the order of the given vector.

Parameters

<i>v</i>	A pointer to the vector to be shuffled.
----------	---

12.19.3.15 `template<typename T> std::vector< T > CX::CX_RandomNumberGenerator::shuffleVector (std::vector< T > v)`

Makes a copy of the given vector, randomizes the order of its elements, and returns the shuffled copy.

Parameters

<i>v</i>	The vector to be operated on.
----------	-------------------------------

Returns

A shuffled copy of v.

12.19.3.16 `double CX_RandomNumberGenerator::uniformDeviate (double lowerBound_closed, double upperBound_open)`

Samples a deviate from a uniform distribution with the range [lowerBound_closed, upperBound_open).

Parameters

<i>lowerBound_ - closed</i>	The lower bound of the distribution. This bound is closed, meaning that you can observe deviates with this value.
<i>upperBound_ - open</i>	The upper bound of the distribution. This bound is open, meaning that you cannot observe deviates with this value.

Returns

The deviate.

12.19.3.17 `std::vector< double > CX_RandomNumberGenerator::uniformDeviates (unsigned int count, double lowerBound_closed, double upperBound_open)`

Samples count deviates from a uniform distribution with the range [*lowerBound_closed*, *upperBound_open*).

Parameters

<i>count</i>	The number of deviates to generate.
<i>lowerBound_ - closed</i>	The lower bound of the distribution. This bound is closed, meaning that you can observe deviates with this value.
<i>upperBound_ - open</i>	The upper bound of the distribution. This bound is open, meaning that you cannot observe deviates with this value.

Returns

A vector of the deviates.

The documentation for this class was generated from the following files:

- CX_RandomNumberGenerator.h
- CX_RandomNumberGenerator.cpp

12.20 CX::CX_SafeDataFrame Class Reference

Inherits [CX::CX_DataFrame](#).

Public Member Functions

- [CX_DataFrameCell](#) **operator()** (std::string column, rowIndex_t row)
- [CX_DataFrameCell](#) **operator()** (rowIndex_t row, std::string column)

Additional Inherited Members

12.20.1 Detailed Description

The documentation for this class was generated from the following files:

- CX_DataFrame.h
- CX_DataFrame.cpp

12.21 CX::CX_Slide_t Struct Reference

```
#include <CX_SlidePresenter.h>
```

Public Types

- enum {
NOT_STARTED, **COPY_TO_BACK_BUFFER_PENDING**, **SWAP_PENDING**, **IN_PROGRESS**,
FINISHED }

Public Attributes

- std::string [slideName](#)
The name of the slide. Set by the user during slide creation.
- ofFbo [framebuffer](#)
A framebuffer containing image data that will be drawn to the screen during this slide's presentation. If drawingFunction points to a user function, framebuffer will not be drawn.
- void(* [drawingFunction](#))(void)
Pointer to a user function that will be called to draw the slide. If this points to a user function, it overrides framebuffer. The drawing function does not need to call ofBackground() or otherwise clear the display before drawing, which allow you to do what is essentially single-buffering using the back buffer as the framebuffer.
- enum CX::CX_Slide_t:: { ... } [slideStatus](#)
- [CX_SlideTimingInfo_t](#) [intended](#)
The intended timing parameters (i.e. what should have happened if there were no presentation errors).
- [CX_SlideTimingInfo_t](#) [actual](#)
The actual timing parameters.
- CX_Micros [copyToBackBufferCompleteTime](#)
The time at which the drawing operations for this slide finished. This is pretty useful to determine if there was an error on the trial (e.g. framebuffer was copied late). If this is greater than actual.startTime, the slide may not have been fully drawn at the time the front and back buffers swapped.

12.21.1 Detailed Description

This struct contains information related to slide presentation using [CX_SlidePresenter](#).

12.21.2 Member Enumeration Documentation

12.21.2.1 anonymous enum

Status of the current slide vis a vis presentation.

12.21.3 Member Data Documentation

12.21.3.1 enum { ... } CX::CX_Slide_t::slideStatus

Status of the current slide vis a vis presentation.

The documentation for this struct was generated from the following file:

- [CX_SlidePresenter.h](#)

12.22 CX::CX_SlidePresenter Class Reference

```
#include <CX_SlidePresenter.h>
```

Public Member Functions

- bool [setup](#) (CX_Display *display)
- bool [setup](#) (const CX_SP_Configuration &config)
- virtual void [update](#) (void)
- void [appendSlide](#) (CX_Slide_t slide)
- void [appendSlideFunction](#) (void(*drawingFunction)(void), CX_Micros slideDuration, std::string slideName="")
- void [beginDrawingNextSlide](#) (CX_Micros slideDuration, std::string slideName="")
- void [endDrawingCurrentSlide](#) (void)
- bool [startSlidePresentation](#) (void)
- void [stopSlidePresentation](#) (void)
- *Stops slide presentation.*
- bool [isPresentingSlides](#) (void)
- *Returns true if slide presentation is in progress, even if the first slide has not yet been presented.*
- void [clearSlides](#) (void)
- std::vector< CX_Slide_t > & [getSlides](#) (void)
- std::vector< CX_Micros > [getActualPresentationDurations](#) (void)
- std::vector< unsigned int > [getActualFrameCounts](#) (void)
- CX_SP_PresentationErrors_t [checkForPresentationErrors](#) (void)

Protected Member Functions

- unsigned int [_calculateFrameCount](#) (CX_Micros duration)
- void [_singleCoreBlockingUpdate](#) (void)
- void [_singleCoreThreadedUpdate](#) (void)
- void [_multiCoreUpdate](#) (void)
- void [_renderCurrentSlide](#) (void)
- void [_waitSyncCheck](#) (void)
- void [_finishPreviousSlide](#) (void)
- void [_handleFinalSlide](#) (void)
- void [_prepareNextSlide](#) (void)

Protected Attributes

- CX_SP_Configuration [_config](#)
- CX_Micros [_hoggingStartTime](#)
- bool [_presentingSlides](#)
- bool [_synchronizing](#)
- unsigned int [_currentSlide](#)
- std::vector< CX_Slide_t > [_slides](#)
- bool [_lastFramebufferActive](#)
- bool [_useFenceSync](#)
- bool [_awaitingFenceSync](#)
- GLsync [_fenceSyncObject](#)

12.22.1 Detailed Description

This class is a very useful abstraction that presents slides (a full display's worth) of visual stimuli for fixed durations. See the [basicChangeDetectionTask.cpp](#), [advancedChangeDetectionTask.cpp](#), and [nBack.cpp](#) examples for the usage of this class.

12.22.2 Member Function Documentation

12.22.2.1 void CX_SlidePresenter::appendSlide (CX_Slide_t *slide*)

Add a fully configured slide to the end of the list of slides. The user code must configure several components of the slide:

- If the framebuffer will be used, the framebuffer must be allocated and drawn to.
- If the drawing function will be used, a valid function pointer must be given. A check is made that either the drawing function is set or the framebuffer is allocated and an error is logged if neither is configured.
- The intended duration must be set.
- The name may be set (optional).

Parameters

<i>slide</i>	The slide to append.
--------------	----------------------

12.22.2.2 void CX_SlidePresenter::appendSlideFunction (void(*)(void) *drawingFunction*, CX_Micros *slideDuration*, std::string *slideName* = " ")

Appends a slide to the slide presenter that will call the given drawing function when it comes time to render the slide to the back buffer. Essentially, the drawing function will be called one frame before the front and back buffers are swapped.

Parameters

<i>drawingFunction</i>	A pointer to a function that will draw the data to the slide. The contents of the back buffer are not clear before this function is called, so the function must clear the background to the desired color.
<i>slideDuration</i>	The amount of time to present the slide for. If this is less than or equal to 0, the slide will be ignored.
<i>slideName</i>	The name of the slide. This can be anything and is purely for the user to use to help identify the slide.

12.22.2.3 void CX_SlidePresenter::beginDrawingNextSlide (CX_Micros *slideDuration*, std::string *slideName* = " ")

Prepares the framebuffer of the next slide for drawing so that any drawing commands given between a call to [beginDrawingNextSlide\(\)](#) and [endDrawingCurrentSlide\(\)](#) will cause stimuli to be drawn to the framebuffer of the slide.

Parameters

<i>slideDuration</i>	The amount of time to present the slide for. If this is less than or equal to 0, the slide will be ignored.
----------------------	---

<i>slideName</i>	The name of the slide. This can be anything and is purely for the user to use to help identify the slide.
------------------	---

12.22.2.4 CX_SP_PresentationErrors_t CX_SlidePresenter::checkForPresentationErrors (void)

Checks the timing data from the last presentation of slides for presentation errors. Currently it checks to see if the intended frame count matches the actual frame count of each slide, which indicates if the duration was correct. It also checks to make sure that the framebuffer was copied to the back buffer before the onset of the slide, which would indicate the potential for vertical tearing.

Returns

A struct with information about the errors that occurred on the last presentation of slides.

Note

If [clearSlides\(\)](#) has been called since the end of the presentation, this does nothing as its data has been cleared. If this function is called during slide presentation, the returned struct will have the presentationErrorsSuccessfullyChecked member set to false and an error will be logged.

12.22.2.5 void CX_SlidePresenter::clearSlides (void)

Clears (deletes) all of the slides contained in the slide presenter and stops presentation, if it was in progress.

12.22.2.6 void CX_SlidePresenter::endDrawingCurrentSlide (void)

Ends drawing to the framebuffer of the slide that is currently being drawn to. See [beginDrawingNextSlide\(\)](#).

12.22.2.7 std::vector< unsigned int > CX_SlidePresenter::getActualFrameCounts (void)

Gets a vector containing the number of frames that each of the slides from the last presentation of slides was presented for. Note that these frame counts may be wrong. If [checkForPresentationErrors\(\)](#) not detect any errors, the frame counts are likely to be right, but there is no guarantee.

Returns

A vector containing the frame counts. The frame count corresponding to the first slide added to the slide presenter will be at index 0.

Note

The frame count of the last slide is meaningless. As far as the slide presenter is concerned, as soon as the last slide is put on the screen, it is done presenting the slides. Because the slide presenter is not responsible for removing the last slide from the screen, it has no idea about the duration of that slide.

12.22.2.8 std::vector< CX_Micros > CX_SlidePresenter::getActualPresentationDurations (void)

Gets a vector containing the durations of the slides from the last presentation of slides. Note that these durations may be wrong. If [checkForPresentationErrors\(\)](#) does not detect any errors, the durations are likely to be right, but there is no guarantee.

Returns

A vector containing the durations. The duration corresponding to the first slide added to the slide presenter will be at index 0.

Note

The duration of the last slide is meaningless. As far as the slide presenter is concerned, as soon as the last slide is put on the screen, it is done presenting the slides. Because the slide presenter is not responsible for removing the last slide from the screen, it has no idea about the duration of that slide.

12.22.2.9 `std::vector< CX_Slide_t > &CX_SlidePresenter::getSlides (void)`

Get a reference to the vector of slides held by the slide presenter. If you modify any of the members of any of the slides, you do so at your own risk. This data is mostly useful in a read-only sort of way (when was that slide presented?).

Returns

A reference to the vector of slides.

12.22.2.10 `bool CX_SlidePresenter::setup (CX_Display * display)`

Set up the slide presenter with the given [CX_Display](#) as the display.

Parameters

<i>display</i>	Pointer to the display to use.
----------------	--------------------------------

Returns

False if there was an error during setup, in which case a message will be logged.

12.22.2.11 `bool CX_SlidePresenter::setup (const CX_SP_Configuration & config)`

Set up the slide presenter using the given configuration.

Parameters

<i>config</i>	The configuration to use.
---------------	---------------------------

Returns

False if there was an error during setup, in which case a message will be logged.

12.22.2.12 `bool CX_SlidePresenter::startSlidePresentation (void)`

Start presenting the slides that are stored in the slide presenter. After this function is called, calls to [update\(\)](#) will advance the state of the slide presentation.

Returns

False if an error was encountered while starting presentation, in which case messages will be logged, true otherwise.

12.22.2.13 `void CX_SlidePresenter::update (void)` `[virtual]`

Updates the state of the slide presenter. If the slide presenter is presenting stimuli, [update\(\)](#) must be called very regularly (at least once per millisecond) in order for the slide presenter to function.

The documentation for this class was generated from the following files:

- `CX_SlidePresenter.h`
- `CX_SlidePresenter.cpp`

12.23 CX::CX_SlideTimingInfo_t Struct Reference

```
#include <CX_SlidePresenter.h>
```

Public Attributes

- uint32_t [startFrame](#)
- uint32_t [frameCount](#)
- CX_Micros [startTime](#)
- CX_Micros [duration](#)

12.23.1 Detailed Description

Contains information about the presentation timing of the slide.

12.23.2 Member Data Documentation

12.23.2.1 CX_Micros CX::CX_SlideTimingInfo_t::duration

Time amount of time the slide was/should have been presented for.

12.23.2.2 uint32_t CX::CX_SlideTimingInfo_t::frameCount

The number of frames the slide was/should be presented for.

12.23.2.3 uint32_t CX::CX_SlideTimingInfo_t::startFrame

The frame on which the slide started/should have started. Can be compared with the value given by `Display.getFrameNumber()`.

12.23.2.4 CX_Micros CX::CX_SlideTimingInfo_t::startTime

The time at which the slide was/should have been started. Can be compared with values from `Clock.getTime()`.

The documentation for this struct was generated from the following file:

- `CX_SlidePresenter.h`

12.24 CX::CX_SoundObject Class Reference

```
#include <CX_SoundObject.h>
```

Public Member Functions

- bool [loadFile](#) (string fileName)
- bool [addSound](#) (string fileName, CX_Micros timeOffset)
- bool [addSound](#) (CX_SoundObject so, CX_Micros timeOffset)
- bool [isReadyToPlay](#) (void)
- bool [isLoadedSuccessfully](#) (void)
- bool [applyGain](#) (float gain, int channel=-1)
- bool [multiplyAmplitudeBy](#) (float amount, int channel=-1)

- float [getPositivePeak](#) (void)
- float [getNegativePeak](#) (void)
- void [setLength](#) (CX_Micros length)
- CX_Micros [getLength](#) (void)
- void [stripLeadingSilence](#) (float tolerance)
- void [addSilence](#) (CX_Micros duration, bool atBeginning)
- void [deleteAmount](#) (CX_Micros duration, bool fromBeginning)
- void [multiplySpeed](#) (float speedMultiplier)
- void [resample](#) (float newSampleRate)
- float [getSampleRate](#) (void)
- bool [setChannelCount](#) (int channels)
- int [getChannelCount](#) (void)
- uint64_t [getTotalSampleCount](#) (void)
- uint64_t [getConcurrentSampleCount](#) (void)
- vector< float > & [getRawDataReference](#) (void)

Public Attributes

- string **name**

12.24.1 Detailed Description

This class is a container for a sound. It can load sound files, manipulate the contents of the sound data, add other sounds to an existing sound at specified offsets.

In order to play a [CX_SoundObject](#), you use a [CX::CX_SoundObjectPlayer](#).

See the [soundObject](#) example for an introduction on how to use this class along with a [CX_SoundObjectPlayer](#).

Note

Nearly all functions of this class should be considered blocking. Many of the operations take quite a while to complete because they are performed on a fairly large vector of sound samples.

12.24.2 Member Function Documentation

12.24.2.1 void CX_SoundObject::addSilence (CX_Micros *duration*, bool *atBeginning*)

Adds the specified amount of silence to the [CX_SoundObject](#) at either the beginning or end.

Parameters

<i>duration</i>	Duration of added silence in microseconds. Dependent on the sample rate of the sound. If the sample rate changes, so does the duration of silence.
<i>atBeginning</i>	If true, silence is added at the beginning of the CX_SoundObject . If false, the silence is added at the end.

12.24.2.2 bool CX_SoundObject::addSound (string *fileName*, CX_Micros *timeOffset*)

Uses [loadFile\(string\)](#) and [addSound\(CX_SoundObject, uint64_t\)](#) to add the given file to the current [CX_SoundObject](#) at the given time offset (in microseconds). See those functions for more information.

Parameters

<i>fileName</i>	Name of the sound file to load.
<i>timeOffset</i>	Time at which to add the new sound.

Returns

Returns true if the new sound was added successfully, false otherwise.

12.24.2.3 bool CX_SoundObject::addSound (CX_SoundObject nso, CX_Micros timeOffset)

Adds the sound data in nso at the time offset. If the sample rates of the sounds differ, nso will be resampled to the sample rate of this CX_SoundObject. If the number of channels of nso does not equal the number of channels of this, an attempt will be made to set the number of channels of nso equal to the number of channels of this CX_SoundObject. The data from nso and this CX_SoundObject are merged by adding the amplitudes of the sounds. The result of the addition is clamped between -1 and 1.

Parameters

<i>nso</i>	A sound object. Must be successfully loaded.
<i>timeOffset</i>	Time at which to add the new sound data in microseconds. Dependent on sample rate.

Returns

True if nso was successfully added to this CX_SoundObject, false otherwise.

12.24.2.4 bool CX_SoundObject::applyGain (float decibels, int channel = -1)

Apply gain to the channel in terms of decibels.

Parameters

<i>decibels</i>	Gain to apply. 0 does nothing. Positive values increase volume, negative values decrease volume. Negative infinity is essentially mute, although see multiplyAmplitudeBy() for a more obvious way to do that same operation.
<i>channel</i>	The channel that the gain should be applied to. If channel is less than 0, the gain is applied to all channels.

12.24.2.5 void CX_SoundObject::deleteAmount (CX_Micros duration, bool fromBeginning)

Deletes the specified amount of sound from the CX_SoundObject from either the beginning or end.

Parameters

<i>duration</i>	Duration of removed sound in microseconds.
<i>fromBeginning</i>	If true, sound is deleted from the beginning of the CX_SoundObject's buffer. If false, the sound is deleted from the end, toward the beginning.

12.24.2.6 CX_Micros CX_SoundObject::getLength (void)

Gets the length, in time, of the sound object.

Returns

The length.

12.24.2.7 float CX_SoundObject::getNegativePeak (void)

Finds the minimum amplitude in the sound object.

Returns

The minimum amplitude.

Note

Amplitudes are between -1 and 1, inclusive.

12.24.2.8 float CX_SoundObject::getPositivePeak (void)

Finds the maximum amplitude in the sound object.

Returns

The maximum amplitude.

Note

Amplitudes are between -1 and 1, inclusive.

12.24.2.9 bool CX_SoundObject::loadFile (string fileName)

Loads a sound file with the given file name into the [CX_SoundObject](#). Any pre-existing data in the sound object is deleted. Some sound file types are supported. Others are not. In the limited testing, mp3 and wav files seem to work well. If the file cannot be loaded, descriptive error messages will be logged.

Parameters

<i>fileName</i>	Name of the sound file to load.
-----------------	---------------------------------

Returns

True if the sound given in the fileName was loaded successfully, false otherwise.

12.24.2.10 bool CX_SoundObject::multiplyAmplitudeBy (float amount, int channel = -1)

Apply gain to the sound. The original value is simply multiplied by the amount and then clamped to be within [-1, 1].

Parameters

<i>amount</i>	The gain that should be applied. A value of 0 mutes the channel. 1 does nothing. 2 doubles the amplitude. -1 inverts the waveform.
<i>channel</i>	The channel that the given multiplier should be applied to. If channel is less than 0, the amplitude multiplier is applied to all channels.

12.24.2.11 void CX_SoundObject::resample (float newSampleRate)

Resamples the audio data stored in the [CX_SoundObject](#) by linear interpolation. Linear interpolation is not the ideal way to resample audio data; some audio fidelity is lost, more so than with other resampling techniques. It is, however, very fast compared to higher-quality methods both in terms of run time and programming time. It has acceptable results, at least when the new sample rate is similar to the old sample rate.

Parameters

<i>newSampleRate</i>	The requested sample rate.
----------------------	----------------------------

12.24.2.12 bool CX_SoundObject::setChannelCount (int *newChannelCount*)

Sets the number of channels of the sound. Depending on the old number of channels (N) and the new number of channels (M), the conversion is performed in different ways. If N == M, nothing happens. If N == 1, each of the M new channels is set equal to the value of the single old channel. If M == 1, the new channel is set equal to the arithmetic average of the N old channels. If (N != 1 && M != 1 && M > N), the first N channels are preserved unchanged and the M - N new channels are set to the arithmetic average of the N old channels. Any other combination of M and N is an error condition.

Parameters

<i>newChannel-Count</i>	The number of channels the CX_SoundObject will have after the conversion.
-------------------------	---

Returns

True if the conversion was successful, false if the attempted conversion is unsupported.

12.24.2.13 void CX_SoundObject::setLength (CX_Micros *length*)

Set the length of the sound to the specified length in microseconds. If the new length is longer than the old length, the new data is zeroed (i.e. set to silence).

12.24.2.14 void CX_SoundObject::stripLeadingSilence (float *tolerance*)

Removes leading "silence" from the sound, where silence is defined by the given tolerance. It is unlikely that the beginning of a sound, even if perceived as silent relative to the rest of the sound, has an amplitude of 0. Therefore, a tolerance of 0 is unlikely to prove useful. Using [getPositivePeak\(\)](#) and/or [getNegativePeak\(\)](#) can help to give a reference amplitude of which some small fraction is perceived as "silent".

Parameters

<i>tolerance</i>	All sound data up to and including the first instance of a sample with an amplitude with an absolute value greater than or equal to tolerance is removed from the sound.
------------------	--

The documentation for this class was generated from the following files:

- CX_SoundObject.h
- CX_SoundObject.cpp

12.25 CX::CX_SoundObjectPlayer Class Reference

```
#include <CX_SoundObjectPlayer.h>
```

Public Member Functions

- bool [setup](#) (CX_SoundObjectPlayerConfiguration_t config)
- bool [play](#) (void)
- bool [startPlayingAt](#) (CX_Micros experimentTime)
- bool [stop](#) (void)

- bool **isPlaying** (void)
- [CX_SoundObjectPlayerConfiguration_t](#) **getConfiguration** (void)
- bool [BLOCKING_setSound](#) ([CX_SoundObject](#) *sound)

12.25.1 Detailed Description

This class is used for playing [CX_SoundObjects](#).

12.25.2 Member Function Documentation

12.25.2.1 bool [CX_SoundObjectPlayer::BLOCKING_setSound](#) ([CX_SoundObject](#) * *sound*)

This function is blocking because the sample rate and number of channels of sound are changed to those of the currently open stream.

Parameters

<i>sound</i>	A pointer to a CX_SoundObject that will be set as the current sound for the CX_SoundObjectPlayer . There are a variety of reasons why the sound could fail to be set as the current sound for the player. If sound was not loaded successfully, this function call fails and an error is logged. If it is not possible to convert the number of channels of sound to the number of channels that the CX_SoundObjectPlayer is configured to use, this function call fails and an error is logged.
--------------	--

Returns

True if sound was successfully set to be the current sound, false otherwise.

12.25.2.2 bool [CX_SoundObjectPlayer::play](#) (void)

Attempts to start playing the current [CX_SoundObject](#) associated with the player.

Returns

True if the sound object associated with the player is [ReadyToPlay\(\)](#), false otherwise.

12.25.2.3 bool [CX_SoundObjectPlayer::setup](#) ([CX_SoundObjectPlayerConfiguration_t](#) *config*)

Configures the [CX_SoundObjectPlayer](#) with the given configuration.

12.25.2.4 bool [CX_SoundObjectPlayer::stop](#) (void)

Stop the currently playing sound object, or, if a playback start was cued, cancel the cued playback.

Returns

Always returns true currently.

The documentation for this class was generated from the following files:

- [CX_SoundObjectPlayer.h](#)
- [CX_SoundObjectPlayer.cpp](#)

12.26 CX::CX_SoundStream Class Reference

```
#include <CX_SoundStream.h>
```

Public Member Functions

- bool **open** (CX_SoundStreamConfiguration_t &config)
- bool **close** (void)
- bool **start** (void)
- bool **stop** (void)
- const
CX_SoundStreamConfiguration_t & **getConfiguration** (void)
- uint64_t **getLastSampleNumber** (void)
- void **setLastSampleNumber** (uint64_t sampleNumber)
- CX_Micros **getStreamLatency** (void)
- bool **hasSwappedSinceLastCheck** (void)
- CX_Micros **getLastSwapTime** (void)
- CX_Micros **estimateNextSwapTime** (void)

Static Public Member Functions

- static std::vector< RtAudio::Api > **getCompiledApis** (void)
- static std::vector< std::string > **convertApisToStrings** (vector< RtAudio::Api > apis)
- static std::string **convertApisToString** (vector< RtAudio::Api > apis, std::string delim="\r\n")
- static std::string **convertApiToString** (RtAudio::Api api)
- static std::vector< std::string > **formatsToStrings** (RtAudioFormat formats)
- static std::string **formatsToString** (RtAudioFormat formats, std::string delim="\r\n")
- static std::vector
< RtAudio::DeviceInfo > **getDeviceList** (RtAudio::Api api)
- static std::string **listDevices** (RtAudio::Api api)

Public Attributes

- ofEvent< CX_SSOutputCallback_t > **outputCallbackEvent**
- ofEvent< CX_SSInputCallback_t > **inputCallbackEvent**

12.26.1 Detailed Description

CX_SoundStream uses RtAudio internally, so you are having problems, you might be able to figure out what is going wrong by checking out the page for it: <http://www.music.mcgill.ca/~gary/rtaudio/index.html>

12.26.2 Member Function Documentation

12.26.2.1 bool CX_SoundStream::close (void)

Closes the sound stream.

Returns

False if an error was encountered while closing the stream, true otherwise.

12.26.2.2 `std::string CX_SoundStream::convertApisToString (vector< RtAudio::Api > apis, std::string delim = "\r\n")`
`[static]`

This helper function converts a vector of `RtAudio::Api` to a string, with the specified delimiter between API names.

Parameters

<i>apis</i>	The vector of RtAudio::Api to convert to string.
<i>delim</i>	The delimiter between elements of the string.

Returns

A string containing the names of the APIs.

12.26.2.3 `std::vector< std::string > CX_SoundStream::convertApisToStrings (vector< RtAudio::Api > apis) [static]`

This helper function converts a vector of RtAudio::Api to a vector of strings, using [convertApiToString\(\)](#) for the conversion.

Parameters

<i>apis</i>	A vector of apis to convert to strings.
-------------	---

Returns

A vector of string names of the apis.

12.26.2.4 `std::string CX_SoundStream::convertApiToString (RtAudio::Api api) [static]`

This helper function converts an RtAudio::Api to a string.

Parameters

<i>api</i>	The api to get a string of.
------------	-----------------------------

Returns

A string of the api name.

12.26.2.5 `CX_Micros CX_SoundStream::estimateNextSwapTime (void)`

Estimate the time at which the next buffer swap will occur.

Returns

The estimated time of next swap. This value can be compared with the result of CX::Instances::Clock.getTime().

12.26.2.6 `std::vector< RtAudio::Api > CX_SoundStream::getCompiledApis (void) [static]`

Get a vector containing a list of all of the APIs for which the RtAudio driver has been compiled to use. If the API you want is not available, you might be able to get it by using a different version of RtAudio.

12.26.2.7 `const CX_SoundStreamConfiguration_t& CX::CX_SoundStream::getConfiguration (void) [inline]`

Gets the configuration that was used on the last call to [open\(\)](#). Because some of the configuration options are only suggestions, this function allows you to check what the actual configuration was.

Returns

A const reference to the configuration struct.

12.26.2.8 `std::vector< RtAudio::DeviceInfo > CX_SoundStream::getDeviceList (RtAudio::Api api) [static]`

For the given api, lists all of the devices on the system that support that api.

Parameters

<i>api</i>	Devices that support this API are scanned.
------------	--

Returns

A machine-readable list of information. See http://www.music.mcgill.ca/~gary/rtaudio/struct-RtAudio_1_1DeviceInfo.html for information about the members of the RtAudio::DeviceInfo struct.

12.26.2.9 CX_Micros CX::CX_SoundStream::getLastSwapTime (void) [inline]

Gets the time at which the last buffer swap occurred.

Returns

This time value can be compared with the result of CX::Instances::Clock.getTime().

12.26.2.10 bool CX_SoundStream::hasSwappedSinceLastCheck (void)

This function checks to see if the audio buffers have been swapped since the last time this function was called.

Returns

True if at least one audio buffer has been swapped out, false if no buffers have been swapped.

12.26.2.11 std::string CX_SoundStream::listDevices (RtAudio::Api *api*) [static]

For the given api, lists all of the devices on the system that support that api. Lots of information about each device is given, like supported sample rates, number of input and output channels, etc.

Parameters

<i>api</i>	Devices that support this API are scanned.
------------	--

Returns

A human-readable formatted string containing the scanned information. Can be printed directly to std::cout or elsewhere.

12.26.2.12 bool CX_SoundStream::open (CX_SoundStreamConfiguration_t & *config*)

Opens the sound stream with the specified configuration. If there was an error during configuration, messages will be logged.

Parameters

<i>config</i>	The configuration settings that are desired. Some of the configuration options are only suggestions, so some of the values that are used may differ from the values that are chosen. In those cases, config is updated based on the actually used settings. You can check the configuration later using getConfiguration() .
---------------	--

Returns

True if configuration appeared to be successful, false otherwise.

Note

Opening the stream does not start it. See [start\(\)](#).

12.26.2.13 bool CX_SoundStream::start (void)

Starts the sound stream. The stream must already be [open\(\)](#).

Returns

False if the stream was not started, true if the stream was started or if it was already running.

12.26.2.14 bool CX_SoundStream::stop (void)

Stop the stream, if is running. If there is an error, a message will be logged.

Returns

False if there was an error, true otherwise.

The documentation for this class was generated from the following files:

- CX_SoundStream.h
- CX_SoundStream.cpp

12.27 CX::CX_SoundStreamConfiguration_t Struct Reference

```
#include <CX_SoundStream.h>
```

Public Attributes

- int [inputChannels](#)
The number of input (e.g. microphone) channels to use. If 0, no input will be used.
- int [outputChannels](#)
The number of output channels to use. Currently only stereo and mono are well-supported.
- int [sampleRate](#)
- unsigned int [bufferSize](#)
- RtAudio::Api [api](#)
- RtAudio::StreamOptions [streamOptions](#)
- int [inputDeviceId](#)
The ID of the desired input device. A value of -1 will cause the system default input device to be used.
- int [outputDeviceId](#)
The ID of the desired output device. A value of -1 will cause the system default output device to be used.

12.27.1 Detailed Description

This struct controls the configuration of the [CX_SoundStream](#).

12.27.2 Member Data Documentation

12.27.2.1 RtAudio::Api CX::CX_SoundStreamConfiguration_t::api

This argument depends on your operating system. Using RtAudio::Api::UNSPECIFIED will attempt to pick a working API from those that are available on your system. The API means the type of software interface to use. For example, on

Windows, you can choose from Windows Direct Sound (DS) and ASIO. ASIO is commonly used with audio recording equipment because it has lower latency whereas DS is more of a consumer-grade interface. The choice of API does not affect how you use this class, but it may affect the performance of sound playback.

See <http://www.music.mcgill.ca/~gary/rtaudio/classRtAudio.html#ac9b6f625da88249d08a8409a9db> for a listing of the APIs. See <http://www.music.mcgill.ca/~gary/rtaudio/classRtAudio.-html#afd0bfa26deae9804e18faff59d0273d9> for the ordering of the APIs.

12.27.2.2 unsigned int CX::CX_SoundStreamConfiguration_t::bufferSize

The size of the audio data buffer to use. A larger buffer size means more latency but also a greater potential for audio glitches (clicks and pops). Buffer size is per channel (i.e. if there are two channels and buffer size is set to 256, the actual buffer size will be 512 samples). Defaults to 4096 samples.

12.27.2.3 int CX::CX_SoundStreamConfiguration_t::sampleRate

The requested sample rate for the input and output channels. If, for the selected device(s), this sample cannot be used, the nearest greater sample rate will be chosen. If there is no greater sample rate, the next lower sample rate will be used.

12.27.2.4 RtAudio::StreamOptions CX::CX_SoundStreamConfiguration_t::streamOptions

See http://www.music.mcgill.ca/~gary/rtaudio/structRtAudio_1_1StreamOptions.-html for more information.

flags must not include RTAUDIO_NONINTERLEAVED: The audio data used by CX is interleaved.

The documentation for this struct was generated from the following file:

- CX_SoundStream.h

12.28 CX::CX_SP_Configuration Struct Reference

```
#include <CX_SlidePresenter.h>
```

Public Types

- enum **SwappingMode** { **SINGLE_CORE_BLOCKING_SWAPS**, **SINGLE_CORE_THREADED_SWAPS**, **MULTI_CORE** }

Public Attributes

- [CX_Display](#) * **display**
A pointer to the display to use.
- std::function< void([CX_FinalSlideFunctionInfo_t](#) &)> **finalSlideCallback**
A pointer to a user function that will be called as soon as the final slide is presented.
- [CX_SP_ErrorMode](#) **errorMode**
- bool **deallocateCompletedSlides**
If true, once a slide has been presented, its framebuffer will be deallocated to conserve memory.
- [CX_Micros](#) **preSwapCPUHoggingDuration**
Only used if swappingMode is a single core mode. The amount of time, before a slide is swapped from the back buffer to the front buffer, that the CPU is put into a spinloop waiting for the buffers to swap.
- enum
CX::CX_SP_Configuration::SwappingMode **swappingMode**

12.28.1 Detailed Description

This struct is used for configuring a [CX_SlidePresenter](#).

The documentation for this struct was generated from the following file:

- [CX_SlidePresenter.h](#)

12.29 CX::CX_SP_PresentationErrors_t Struct Reference

```
#include <CX_SlidePresenter.h>
```

Public Member Functions

- unsigned int [totalErrors](#) (void)
Sums up all of the different types of errors that are measured.

Public Attributes

- bool [presentationErrorsSuccessfullyChecked](#)
True if presentation errors were successfully checked for. This does not mean that there were no presentation errors, but that there were no presentation error checking errors.
- unsigned int [incorrectFrameCounts](#)
The number of slides for which the actual and intended frame counts did not match.
- unsigned int [lateCopiesToBackBuffer](#)
The number of slides for which the time at which the slide finished being copied to the back buffer was after the actual start time of the slide.

12.29.1 Detailed Description

Contains information about errors that were detected during slide presentation. See [CX_SlidePresenter::checkForPresentationErrors\(\)](#).

The documentation for this struct was generated from the following file:

- [CX_SlidePresenter.h](#)

12.30 CX::Util::CX_TrialController Class Reference

```
#include <CX_TrialController.h>
```

Public Member Functions

- int [update](#) (void)
- void [start](#) (void)
"Arm" the trial controller. Before this is called, [update\(\)](#) will do nothing.
- void [stop](#) (void)
"Disarm" the trial controller. After this is called, [update\(\)](#) will do nothing.
- bool [isActive](#) (void)

Check to see if the trial controller is active. See [start\(\)](#) and [stop\(\)](#).

- void [appendFunction](#) (std::function< int(void)> userFunction)
- void [reset](#) (void)
- bool [setCurrentFunction](#) (int currentFunction)
- int [getCurrentFunction](#) (void)

Get the index of the current function (i.e. the function that will be called the next time [update\(\)](#) is called).

- unsigned int [getFunctionCount](#) (void)

Get the number of user functions stored by this trial controller.

12.30.1 Detailed Description

This class is used to help with the fact that most psychology experiments are by nature more or less linear, but that [CX](#) works better if code does not block (see [Blocking Code](#)).

The way this works is that segments of user code, each representing one part of a trial, are put into functions. Those functions are added to the trial controller with [appendFunction\(\)](#), which puts the function at the end of the list of functions. User functions take no arguments and return an int.

When you want to use the trial controller, call [start\(\)](#) and it will be "armed". When it is armed and [update\(\)](#) is called, it will call the current function in the list of user functions. If the user function is done with whatever it needs to do, it should return 1. This will cause the trial controller to move on to the next user function. If the user function is not done with its task, it should return 0. If it returns 0, it will be called again the next time [update\(\)](#) is called.

12.30.2 Member Function Documentation

12.30.2.1 void CX_TrialController::appendFunction (std::function< int(void)> userFunction)

Adds a user function to the end of the list of functions to be called by the trial controller.

Parameters

<i>userFunction</i>	Typically a pointer to a function that takes no arguments and returns an int. Because it is a std::function, it can also be a lamda.
---------------------	--

12.30.2.2 void CX_TrialController::reset (void)

Clear the user functions and otherwise reset to default state.

Note

This function stops the trial controller ([isActive\(\)](#) will return false).

12.30.2.3 bool CX_TrialController::setCurrentFunction (int currentFunction)

Sets the current user function by index, which allows you to skip over functions or go back to a previous function.

Parameters

<i>currentFunction</i>	The new current function index. If this is out of range, an error will be logged and the function will return false.
------------------------	--

Returns

False if the index was out of range, true otherwise.

Note

If this is called from within a user function that has been called from this instance of the [CX_TrialController](#), that function should return 0. If it does not, [setCurrentFunction\(\)](#) will set the function index and then that index will be incremented after the user function completes. However, if 0 is returned, the function index is not incremented.

12.30.2.4 int CX_TrialController::update (void)

Updates the trial controller state. Each time this function is called, the user function at the current function index is called. If that function returns a nonzero value, the trial controller will increment the current function index and that function will be called the next time [update\(\)](#) is called. If the current function index is incremented past the end of the list of functions, it will wrap around to the beginning of the list.

Returns

The value returned by the user function that was called.

Note

This function should probably be called every time `updateExperiment()` is called, although there are other use cases.

If not [isActive\(\)](#), this function does nothing and returns 0.

The documentation for this class was generated from the following files:

- CX_TrialController.h
- CX_TrialController.cpp

13 File Documentation

13.1 advancedChangeDetectionTask.cpp File Reference

```
#include "CX_EntryPoint.h"
```

Functions

- int **drawStimuli** (void)
- int **presentStimuli** (void)
- int **getResponse** (void)
- void **generateTrials** (int trialCount)
- void **updateExperiment** (void)
- void **drawFixation** (void)
- void **drawBlank** (void)
- void **drawSampleArray** (void)
- void **drawTestArray** (void)
- ofColor **backgroundColor** (50)
- void [runExperiment](#) (void)

Variables

- [CX_TrialController](#) **trialController**
- [CX_SlidePresenter](#) **SlidePresenter**
- [CX_SafeDataFrame](#) **trialDf**
- int **trialIndex** = 0
- int **circleRadius** = 30

13.1.1 Detailed Description

This example is a more advanced version of the change detection task presented in the `basicChangeDetectionTask` example. It is not "advanced" because it is more complex, but because it uses more features of [CX](#). It actually ends up being more simple because of how it uses features of [CX](#).

Items that are commented are new, although not all new stuff will necessarily be commented. The two main features that are demonstrated are `CX_SafeDataFrame` and `CX_TrialController`.

13.2 `basicChangeDetectionTask.cpp` File Reference

```
#include "CX_EntryPoint.h"
```

Functions

- void **updateExperiment** (void)
- vector< TrialData_t > **generateTrials** (int trialCount)
- void **outputData** (void)
- void **drawFixation** (void)
- void **drawBlank** (void)
- void **drawSampleArray** (const TrialData_t &tr)
- void **drawTestArray** (const TrialData_t &tr)
- ofColor **backgroundColor** (50)
- void [runExperiment](#) (void)

Variables

- [CX_SlidePresenter](#) **SlidePresenter**
- vector< TrialData_t > **trials**
- int **trialIndex** = 0
- int **circleRadius** = 30
- string **trialPhase** = "drawStimuli"

13.2.1 Detailed Description

This example shows how to do a simple change-detection experiment using [CX](#). The stimuli are colored circles which are presented in a 3X3 matrix.

Press the S key to indicate that you think the test array is the same or the D key to indicate that you think the test array is different.

13.3 nBack.cpp File Reference

```
#include "CX_EntryPoint.h"
```

Functions

- ofColor **backgroundColor** (50)
- ofColor **textColor** (255)
- void **lastSlideFunction** ([CX_FinalSlideFunctionInfo_t](#) &info)
- void **drawStimulusForTrial** (unsigned int trial, bool showInstructions)
- void **generateTrials** (int numberOfTrials)
- void [runExperiment](#) (void)

Variables

- [CX_DataFrame](#) **df**
- [CX_DataFrame::rowIndex_t](#) **trialNumber** = 0
- int **trialCount** = 40
- int **lastSlideIndex** = 0
- int **nBack** = 2
- ofTrueTypeFont **letterFont**
- ofTrueTypeFont **instructionFont**
- char **targetKey** = 'f'
- char **nonTargetKey** = 'j'
- [CX_Millis](#) **stimulusPresentationDuration** = 1000.0
- [CX_Millis](#) **interStimulusInterval** = 1000.0
- [CX_SlidePresenter](#) **SlidePresenter**

13.3.1 Detailed Description

This example shows how to implement an N-Back task using an advanced feature of the [CX_SlidePresenter](#) (SP). There is a feature of the SP that allows you to give it a pointer to a function that will be called every time the SP has just presented the final slide that it currently has. In your function, you can add more slides to the SP, which will allow it to continue presenting slides. If you don't add any more slides, slide presentation will stop with the currently presented slide.

For this N-Back task, the presentation of stimuli will follow the pattern stimulus-blank-stimulus-blank etc. The idea is that you will load up the SP with the first few stimuli and blanks. The SP will be started and will present the first few stimuli. When it runs out of stimuli, the last slide user function will be called. In this function, we will check for any responses that have been made since the last time the function was called and draw the next stimulus-blank pair. See the definition of [lastSlideFunction](#) and [setupExperiment](#) for the implementation of these ideas.

Index

AXIS_POSITION_CHANGE
 CX::CX_JoystickEvent_t, 52
addColumn
 CX::CX_DataFrame, 31
addSilence
 CX::CX_SoundObject, 74
addSound
 CX::CX_SoundObject, 74, 75
advancedChangeDetectionTask.cpp, 87
api
 CX::CX_SoundStreamConfiguration_t, 83
appendFunction
 CX::Util::CX_TrialController, 86
appendRow
 CX::CX_DataFrame, 32
appendSlide
 CX::CX_SlidePresenter, 70
appendSlideFunction
 CX::CX_SlidePresenter, 70
applyGain
 CX::CX_SoundObject, 75
arrayToVector
 CX::Util, 22
at
 CX::CX_DataFrame, 32
availableEvents
 CX::CX_Joystick, 50
 CX::CX_Keyboard, 53
 CX::CX_Mouse, 58

BUTTON_PRESS
 CX::CX_JoystickEvent_t, 52
BUTTON_RELEASE
 CX::CX_JoystickEvent_t, 52
BLOCKING_estimateFramePeriod
 CX::CX_Display, 44
BLOCKING_setAutoSwapping
 CX::CX_Display, 44
BLOCKING_setSound
 CX::CX_SoundObjectPlayer, 78
BLOCKING_swapFrontAndBackBuffers
 CX::CX_Display, 44
BLOCKING_waitForOpenGL
 CX::CX_Display, 44
basicChangeDetectionTask.cpp, 88
beginDrawingNextSlide
 CX::CX_SlidePresenter, 70
beginDrawingToBackBuffer
 CX::CX_Display, 44
binomialDeviate
 CX::CX_RandomNumberGenerator, 61

bufferSize
 CX::CX_SoundStreamConfiguration_t, 84

CX
 FIX_TIMING_FROM_FIRST_SLIDE, 18
 PROPAGATE_DELAYS, 18
CX::CX_JoystickEvent_t
 AXIS_POSITION_CHANGE, 52
 BUTTON_PRESS, 52
 BUTTON_RELEASE, 52
CX::CX_KeyEvent_t
 PRESSED, 54
 RELEASED, 54
 REPEAT, 54
CX::CX_MouseEvent_t
 DRAGGED, 60
 MOVED, 60
 PRESSED, 60
 RELEASED, 60
 SCROLLED, 60
CX, 17
 CX_SP_ErrorMode, 18
CX::CX_Clock, 26
 getDateTimeString, 26
 getExperimentStartDateTimeString, 27
 getExperimentStartTime, 27
 getSystemTime, 27
 getTime, 27
CX::CX_DataFrame, 30
 addColumn, 31
 appendRow, 32
 at, 32
 clear, 32
 columnNames, 32
 copyColumn, 32
 copyColumns, 32
 copyRows, 34
 deleteColumn, 34
 deleteRow, 34
 operator(), 34
 operator=, 35
 print, 35
 printToFile, 36
 readFromFile, 36
 reorderRows, 36
 setRowCount, 37
 shuffleRows, 37
CX::CX_DataFrameCell, 37
 CX_DataFrameCell, 39
 copyCellTo, 39
 getStoredType, 39

- operator=, 39
- store, 39
- storeVector, 39
- to, 39, 40
- toVector, 40
- CX::CX_DataFrameColumn, 40
- CX::CX_DataFrameRow, 41
- CX::CX_Display, 43
 - BLOCKING_estimateFramePeriod, 44
 - BLOCKING_setAutoSwapping, 44
 - BLOCKING_swapFrontAndBackBuffers, 44
 - BLOCKING_waitForOpenGL, 44
 - beginDrawingToBackBuffer, 44
 - drawFboToBackBuffer, 44
 - endDrawingToBackBuffer, 46
 - estimateNextSwapTime, 46
 - getCenterOfDisplay, 46
 - getFrameNumber, 46
 - getFramePeriod, 46
 - getLastSwapTime, 46
 - getResolution, 46
 - hasSwappedSinceLastCheck, 46
 - isAutomaticallySwapping, 47
 - setFullScreen, 47
 - setWindowResolution, 47
 - setWindowTitle, 47
 - setup, 47
 - swapFrontAndBackBuffers, 47
- CX::CX_FinalSlideFunctionInfo_t, 47
- CX::CX_InputManager, 48
 - pollEvents, 48
 - setup, 48
- CX::CX_Joystick, 50
 - availableEvents, 50
 - clearEvents, 50
 - getAxisPositions, 50
 - getButtonStates, 50
 - getJoystickName, 51
 - getNextEvent, 51
 - pollEvents, 51
 - setup, 51
- CX::CX_JoystickEvent_t, 51
 - JoystickEventType, 52
- CX::CX_KeyEvent_t, 53
 - key, 54
 - KeyboardEventType, 54
- CX::CX_Keyboard, 52
 - availableEvents, 53
 - clearEvents, 53
 - getNextEvent, 53
- CX::CX_Logger, 55
 - captureOFLogMessages, 56
 - error, 56
 - fatalError, 56
 - flush, 56
 - level, 56
 - levelForAllModules, 56
 - levelForConsole, 56
 - levelForFile, 56
 - log, 57
 - notice, 57
 - setMessageFlushCallback, 57
 - timestamps, 57
 - verbose, 57
 - warning, 57
- CX::CX_Mouse, 58
 - availableEvents, 58
 - clearEvents, 58
 - getCursorPosition, 58
 - getNextEvent, 58
 - setCursorPosition, 58
 - showCursor, 59
- CX::CX_MouseEvent_t, 59
 - MouseEventType, 60
- CX::CX_RandomNumberGenerator, 60
 - binomialDeviates, 61
 - CX_RandomNumberGenerator, 61
 - getMaximumRandomInt, 61
 - getMinimumRandomInt, 62
 - getSeed, 62
 - normalDeviates, 62
 - randomExclusive, 62
 - randomInt, 64
 - sample, 64, 65
 - setSeed, 65
 - shuffleVector, 65
 - uniformDeviate, 65
 - uniformDeviates, 67
- CX::CX_SP_Configuration, 84
- CX::CX_SP_PresentationErrors_t, 85
- CX::CX_SafeDataFrame, 67
- CX::CX_Slide_t, 68
 - slideStatus, 68
- CX::CX_SlidePresenter, 69
 - appendSlide, 70
 - appendSlideFunction, 70
 - beginDrawingNextSlide, 70
 - checkForPresentationErrors, 71
 - clearSlides, 71
 - endDrawingCurrentSlide, 71
 - getActualFrameCounts, 71
 - getActualPresentationDurations, 71
 - getSlides, 72
 - setup, 72
 - startSlidePresentation, 72
 - update, 72
- CX::CX_SlideTimingInfo_t, 73
 - duration, 73

- frameCount, 73
- startFrame, 73
- startTime, 73
- CX::CX_SoundObject, 73
 - addSilence, 74
 - addSound, 74, 75
 - applyGain, 75
 - deleteAmount, 75
 - getLength, 75
 - getNegativePeak, 75
 - getPositivePeak, 76
 - loadFile, 76
 - multiplyAmplitudeBy, 76
 - resample, 76
 - setChannelCount, 77
 - setLength, 77
 - stripLeadingSilence, 77
- CX::CX_SoundObjectPlayer, 77
 - BLOCKING_setSound, 78
 - play, 78
 - setup, 78
 - stop, 78
- CX::CX_SoundStream, 79
 - close, 79
 - convertApiToString, 81
 - convertApisToString, 79
 - convertApisToStrings, 81
 - estimateNextSwapTime, 81
 - getCompiledApis, 81
 - getConfiguration, 81
 - getDeviceList, 81
 - getLastSwapTime, 82
 - hasSwappedSinceLastCheck, 82
 - listDevices, 82
 - open, 82
 - start, 82
 - stop, 83
- CX::CX_SoundStreamConfiguration_t, 83
 - api, 83
 - bufferSize, 84
 - sampleRate, 84
 - streamOptions, 84
- CX::Draw, 18
 - centeredString, 19
 - squircleToPath, 19
 - star, 19
 - starToPath, 20
- CX::Instances, 20
- CX::Private, 20
- CX::Util, 21
 - arrayToVector, 22
 - checkOFVersion, 22
 - intVector, 22
 - pixelsToDegrees, 22
 - repeat, 22, 24
 - round, 24
 - sequence, 25
 - sequenceAlong, 25
 - sequenceSteps, 25
 - writeToFile, 25
- CX::Util::CX_CoordinateConverter, 27
 - CX_CoordinateConverter, 28
 - operator(), 28, 29
 - setAxisInversion, 29
 - setOrigin, 29
 - setUnitConverter, 29
- CX::Util::CX_DegreeToPixelConverter, 41
 - CX_DegreeToPixelConverter, 42
 - operator(), 43
- CX::Util::CX_LengthToPixelConverter, 54
 - CX_LengthToPixelConverter, 54
 - operator(), 55
- CX::Util::CX_TrialController, 85
 - appendFunction, 86
 - reset, 86
 - setCurrentFunction, 86
 - update, 87
- CX_CoordinateConverter
 - CX::Util::CX_CoordinateConverter, 28
- CX_DataFrameCell
 - CX::CX_DataFrameCell, 39
- CX_DegreeToPixelConverter
 - CX::Util::CX_DegreeToPixelConverter, 42
- CX_LengthToPixelConverter
 - CX::Util::CX_LengthToPixelConverter, 54
- CX_LogLevel
 - Error Logging, 10
- CX_RandomNumberGenerator
 - CX::CX_RandomNumberGenerator, 61
- CX_SP_ErrorMode
 - CX, 18
- captureOFLogMessages
 - CX::CX_Logger, 56
- centeredString
 - CX::Draw, 19
- checkForPresentationErrors
 - CX::CX_SlidePresenter, 71
- checkOFVersion
 - CX::Util, 22
- clear
 - CX::CX_DataFrame, 32
- clearEvents
 - CX::CX_Joystick, 50
 - CX::CX_Keyboard, 53
 - CX::CX_Mouse, 58
- clearSlides
 - CX::CX_SlidePresenter, 71
- Clock

- Timing, [14](#)
- close
 - CX::CX_SoundStream, [79](#)
- columnNames
 - CX::CX_DataFrame, [32](#)
- convertApiToString
 - CX::CX_SoundStream, [81](#)
- convertApisToString
 - CX::CX_SoundStream, [79](#)
- convertApisToStrings
 - CX::CX_SoundStream, [81](#)
- copyCellTo
 - CX::CX_DataFrameCell, [39](#)
- copyColumn
 - CX::CX_DataFrame, [32](#)
- copyColumns
 - CX::CX_DataFrame, [32](#)
- copyRows
 - CX::CX_DataFrame, [34](#)
- DRAGGED
 - CX::CX_MouseEvent_t, [60](#)
- Data, [8](#)
- deleteAmount
 - CX::CX_SoundObject, [75](#)
- deleteColumn
 - CX::CX_DataFrame, [34](#)
- deleteRow
 - CX::CX_DataFrame, [34](#)
- Display
 - Entry Point, [9](#)
- drawFboToBackBuffer
 - CX::CX_Display, [44](#)
- duration
 - CX::CX_SlideTimingInfo_t, [73](#)
- endDrawingCurrentSlide
 - CX::CX_SlidePresenter, [71](#)
- endDrawingToBackBuffer
 - CX::CX_Display, [46](#)
- Entry Point, [9](#)
 - Display, [9](#)
 - Input, [9](#)
 - Log, [9](#)
 - RNG, [9](#)
 - runExperiment, [9](#)
- error
 - CX::CX_Logger, [56](#)
- Error Logging, [10](#)
 - CX_LogLevel, [10](#)
- estimateNextSwapTime
 - CX::CX_Display, [46](#)
 - CX::CX_SoundStream, [81](#)
- FIX_TIMING_FROM_FIRST_SLIDE
 - CX, [18](#)
- fatalError
 - CX::CX_Logger, [56](#)
- flush
 - CX::CX_Logger, [56](#)
- frameCount
 - CX::CX_SlideTimingInfo_t, [73](#)
- getActualFrameCounts
 - CX::CX_SlidePresenter, [71](#)
- getActualPresentationDurations
 - CX::CX_SlidePresenter, [71](#)
- getAxisPositions
 - CX::CX_Joystick, [50](#)
- getButtonStates
 - CX::CX_Joystick, [50](#)
- getCenterOfDisplay
 - CX::CX_Display, [46](#)
- getCompiledApis
 - CX::CX_SoundStream, [81](#)
- getConfiguration
 - CX::CX_SoundStream, [81](#)
- getCursorPosition
 - CX::CX_Mouse, [58](#)
- getDateTimeString
 - CX::CX_Clock, [26](#)
- getDeviceList
 - CX::CX_SoundStream, [81](#)
- getExperimentStartDateTimeString
 - CX::CX_Clock, [27](#)
- getExperimentStartTime
 - CX::CX_Clock, [27](#)
- getFrameNumber
 - CX::CX_Display, [46](#)
- getFramePeriod
 - CX::CX_Display, [46](#)
- getJoystickName
 - CX::CX_Joystick, [51](#)
- getLastSwapTime
 - CX::CX_Display, [46](#)
 - CX::CX_SoundStream, [82](#)
- getLength
 - CX::CX_SoundObject, [75](#)
- getMaximumRandomInt
 - CX::CX_RandomNumberGenerator, [61](#)
- getMinimumRandomInt
 - CX::CX_RandomNumberGenerator, [62](#)
- getNegativePeak
 - CX::CX_SoundObject, [75](#)
- getNextEvent
 - CX::CX_Joystick, [51](#)
 - CX::CX_Keyboard, [53](#)
 - CX::CX_Mouse, [58](#)
- getPositivePeak

- CX::CX_SoundObject, 76
- getResolution
 - CX::CX_Display, 46
- getSeed
 - CX::CX_RandomNumberGenerator, 62
- getSlides
 - CX::CX_SlidePresenter, 72
- getStoredType
 - CX::CX_DataFrameCell, 39
- getSystemTime
 - CX::CX_Clock, 27
- getTime
 - CX::CX_Clock, 27
- hasSwappedSinceLastCheck
 - CX::CX_Display, 46
 - CX::CX_SoundStream, 82
- Input
 - Entry Point, 9
- Input Devices, 11
- intVector
 - CX::Util, 22
- isAutomaticallySwapping
 - CX::CX_Display, 47
- JoystickEventType
 - CX::CX_JoystickEvent_t, 52
- key
 - CX::CX_KeyEvent_t, 54
- KeyboardEventType
 - CX::CX_KeyEvent_t, 54
- level
 - CX::CX_Logger, 56
- levelForAllModules
 - CX::CX_Logger, 56
- levelForConsole
 - CX::CX_Logger, 56
- levelForFile
 - CX::CX_Logger, 56
- listDevices
 - CX::CX_SoundStream, 82
- loadFile
 - CX::CX_SoundObject, 76
- Log
 - Entry Point, 9
- log
 - CX::CX_Logger, 57
- MOVED
 - CX::CX_MouseEvent_t, 60
- MouseEventType
 - CX::CX_MouseEvent_t, 60
- multiplyAmplitudeBy
 - CX::CX_SoundObject, 76
- nBack.cpp, 89
- normalDeviates
 - CX::CX_RandomNumberGenerator, 62
- notice
 - CX::CX_Logger, 57
- open
 - CX::CX_SoundStream, 82
- operator()
 - CX::CX_DataFrame, 34
 - CX::Util::CX_CoordinateConverter, 28, 29
 - CX::Util::CX_DegreeToPixelConverter, 43
 - CX::Util::CX_LengthToPixelConverter, 55
- operator=
 - CX::CX_DataFrame, 35
 - CX::CX_DataFrameCell, 39
- PRESSED
 - CX::CX_KeyEvent_t, 54
 - CX::CX_MouseEvent_t, 60
- PROPAGATE_DELAYS
 - CX, 18
- pixelsToDegrees
 - CX::Util, 22
- play
 - CX::CX_SoundObjectPlayer, 78
- pollEvents
 - CX::CX_InputManager, 48
 - CX::CX_Joystick, 51
- print
 - CX::CX_DataFrame, 35
- printToFile
 - CX::CX_DataFrame, 36
- RELEASED
 - CX::CX_KeyEvent_t, 54
 - CX::CX_MouseEvent_t, 60
- REPEAT
 - CX::CX_KeyEvent_t, 54
- RNG
 - Entry Point, 9
- randomExclusive
 - CX::CX_RandomNumberGenerator, 62
- randomInt
 - CX::CX_RandomNumberGenerator, 64
- Randomization, 12
- readFromFile
 - CX::CX_DataFrame, 36
- reorderRows
 - CX::CX_DataFrame, 36
- repeat
 - CX::Util, 22, 24

- resample
 - CX::CX_SoundObject, 76
- reset
 - CX::Util::CX_TrialController, 86
- round
 - CX::Util, 24
- runExperiment
 - Entry Point, 9
- SCROLLED
 - CX::CX_MouseEvent_t, 60
- sample
 - CX::CX_RandomNumberGenerator, 64, 65
- sampleRate
 - CX::CX_SoundStreamConfiguration_t, 84
- sequence
 - CX::Util, 25
- sequenceAlong
 - CX::Util, 25
- sequenceSteps
 - CX::Util, 25
- setAxisInversion
 - CX::Util::CX_CoordinateConverter, 29
- setChannelCount
 - CX::CX_SoundObject, 77
- setCurrentFunction
 - CX::Util::CX_TrialController, 86
- setCursorPosition
 - CX::CX_Mouse, 58
- setFullScreen
 - CX::CX_Display, 47
- setLength
 - CX::CX_SoundObject, 77
- setMessageFlushCallback
 - CX::CX_Logger, 57
- setOrigin
 - CX::Util::CX_CoordinateConverter, 29
- setRowCount
 - CX::CX_DataFrame, 37
- setSeed
 - CX::CX_RandomNumberGenerator, 65
- setUnitConverter
 - CX::Util::CX_CoordinateConverter, 29
- setWindowResolution
 - CX::CX_Display, 47
- setWindowTitle
 - CX::CX_Display, 47
- setup
 - CX::CX_Display, 47
 - CX::CX_InputManager, 48
 - CX::CX_Joystick, 51
 - CX::CX_SlidePresenter, 72
 - CX::CX_SoundObjectPlayer, 78
- showCursor
 - CX::CX_Mouse, 59
- shuffleRows
 - CX::CX_DataFrame, 37
- shuffleVector
 - CX::CX_RandomNumberGenerator, 65
- slideStatus
 - CX::CX_Slide_t, 68
- Sound, 13
- squircleToPath
 - CX::Draw, 19
- star
 - CX::Draw, 19
- starToPath
 - CX::Draw, 20
- start
 - CX::CX_SoundStream, 82
- startFrame
 - CX::CX_SlideTimingInfo_t, 73
- startSlidePresentation
 - CX::CX_SlidePresenter, 72
- startTime
 - CX::CX_SlideTimingInfo_t, 73
- stop
 - CX::CX_SoundObjectPlayer, 78
 - CX::CX_SoundStream, 83
- store
 - CX::CX_DataFrameCell, 39
- storeVector
 - CX::CX_DataFrameCell, 39
- streamOptions
 - CX::CX_SoundStreamConfiguration_t, 84
- stripLeadingSilence
 - CX::CX_SoundObject, 77
- swapFrontAndBackBuffers
 - CX::CX_Display, 47
- timestamps
 - CX::CX_Logger, 57
- Timing, 14
 - Clock, 14
- to
 - CX::CX_DataFrameCell, 39, 40
- toVector
 - CX::CX_DataFrameCell, 40
- uniformDeviate
 - CX::CX_RandomNumberGenerator, 65
- uniformDeviates
 - CX::CX_RandomNumberGenerator, 67
- update
 - CX::CX_SlidePresenter, 72
 - CX::Util::CX_TrialController, 87
- Utility, 15
- verbose

CX::CX_Logger, [57](#)
Video, [16](#)

warning
 CX::CX_Logger, [57](#)
writeToFile
 CX::Util, [25](#)